

Davin

Orthopedic Business
Worldwide Headquarters

NO. 123 Guangyuan Road (Building 21-2)
Zouqu Town, Zhonglou District
Changzhou City

CONTACT NUMBER**+86 18932382609**

The product technique shown is for illustrative purposes only. The technique(s) actually employed in each case will always depend upon the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient.

Please see the package insert for the complete list of indications, warnings, precautions, and other important medical information.

DAVIN Trauma PRODUCT MANUAL

People-oriented technology



For more information visit

<https://dafinki.en.alibaba.com>

www.Davin Med.com

COMPANY PROFILE



Davin Medical is a professional manufacturer specializing in the production of Class III orthopedic implants. It has 6 years of R&D and manufacturing experience. A professional overseas trade team provides you with product solutions and market growth services.

We offer a wide range of implants, such as plates, screws, intramedullary nails, spinal products, and orthopedic instruments. Our products have obtained CE and ISO 13485 certifications.

We have more than 100 employees, including over 30 R&D personnel and more than 20 domestic and international sales staff.



Category	Annual Capability
Spinal Screws	100k EA
Bone Screws	1000k EA
Intramedullary Nails	120k EA

纵切中心
CNC Turning Center
CITIZEN from Japan



Category	Annual Capability
Bone Plates	150k EA
External Fixators	20k EA

加工中心
CNC Machining Center
HAAS from USA




FAQ

Q1: What is your payment method?

Answer: We accept bank transfers, Western Union, PayPal, and Alibaba.

Q2: What is your delivery time?

Answer: We have a large inventory and can ship as soon as possible within a week.

Q3: What is the delivery method?

Answer: We ship orders through DHL, UPS, FEDEX, TNT, EMS and other methods.

Q4: What is your product warranty?

Answer: Our product comes with a one-year warranty.

Q5: Do you provide free samples?

Answer: Yes, we can provide free samples, but you must pay for shipping.

指骨锁定板 Phalange Locking Plates

指骨直型锁定板 Straight Phalange Locking Plates

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3092N040	4孔 4H	23	纯钛 Pure Titanium
	KZ011-3092N060	6孔 6H	36	
	备注: 选用1.5锁定螺钉, 1.5皮质骨螺钉 Remark:Select 1.5 locking screws,1.5 cortical screws			

指骨重建型锁定板 Phalange Reconstruction Locking Plate



Product No.	Specification	Length	Material
KZ011-3102N040	4孔 4H	19	纯钛 Pure Titanium
KZ011-3102N060	6孔 6H	29	
KZ011-3102N070	7孔 7H	34	
KZ011-3102N080	8孔 8H	39	
KZ011-3102N090	9孔 9H	44	
KZ011-3102N100	10孔 10H	49	
KZ011-3102N110	11孔 11H	54	
KZ011-3102N120	12孔 12H	59	
备注：选用1.5锁定螺钉 Remark:Select 1.5 locking screws			

指骨支撑型锁定板 Phalange Support Locking Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3142L080	8孔 8H	23	纯钛 Pure Titanium
	KZ011-3142R080	8孔 8H	23	
备注: 选用1.5锁定螺钉 Remark:Select 1.5 locking screws				

指骨重建T型锁定板-头三孔（翘） Phalange Reconstruct the T-type locking plate - head three holes (tilted)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3112N040	4孔 4H	26	纯钛 Pure Titanium
	KZ011-3112N050	5孔 5H	31	
	KZ011-3112N060	6孔 6H	36	
	KZ011-3112N070	7孔 7H	41	
	KZ011-3112N080	8孔 8H	46	
备注: 选用1.5锁定螺钉 Remark:Select 1.5 locking screws				

指骨锁定板 Phalange Locking Plates

指骨重建T型锁定板-III Phalange Reconstruct the T-type locking plate-III

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3152N070	7孔 7H	39	纯钛 Pure Titanium
	KZ011-3152N080	8孔 8H	44.5	
	备注: 选用1.5锁定螺钉, 1.5皮质骨螺钉 Remark:Select 1.5 locking screws,1.5 cortical screws			

指骨Y型锁定板 Phalange Y-type locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3042N030	3孔 3H	27	纯钛 Pure Titanium
KZ011-3042N040	4孔 4H	34	
KZ011-3042N070	7孔 7H	40	
KZ011-3042N080	8孔 8H	47	
KZ011-3042N090	9孔 9H	54	
KZ011-3042N100	10孔 10H	61	
备注: 选用2.0锁定螺钉, 2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨重建T型锁定板- II Phalange Reconstruct the T-type locking plate- II


	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3022N030	3孔 3H	25	纯钛 Pure Titanium
	KZ011-3022N040	4孔 4H	32.1	
	KZ011-3022N050	5孔 5H	38.9	
	KZ011-3022N060	6孔 6H	47	
	KZ011-3022N070	7孔 7H	52	
	备注: 选用2.0锁定螺钉, 2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指髌骨锁定板 Phalange Condylar locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3162N040	4孔 4H	36	纯钛 Pure Titanium
	KZ011-3162N050	5孔 5H	43	
	KZ011-3162N060	6孔 6H	47.2	
	KZ011-3162N070	7孔 7H	54	
	备注: 选用2.0锁定螺钉, 2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨锁定板 Phalange Locking Plates

指骨重建T型锁定板- I Phalange Reconstruct the T-type locking plate- I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3032N030	3孔 3H	25	纯钛 Pure Titanium
KZ011-3032N040	4孔 4H	32	
KZ011-3032N050	5孔 5H	39	
KZ011-3032N060	6孔 6H	47	
KZ011-3032N070	7孔 7H	52	
备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨曲型锁定板 Phalange Curved locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3132N070	7孔 7H	46	纯钛 Pure Titanium
	备注：选用2.0锁定螺钉 Remark:Select 2.0 locking screws			

掌髌骨锁定板- I Metacarpal condyle locking plate- I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3292N040	4孔 4H	27	纯钛 Pure Titanium
	KZ011-3292N050	5孔 5H	32	
	备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

掌髌骨锁定板- II Metacarpal condyle locking plate- II

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3052N040	4孔 4H	34	纯钛 Pure Titanium
	KZ011-3052N050	5孔 5H	40	
	备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨T型锁定板 Phalange T-type locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3082N050	5孔 5H	32	纯钛 Pure Titanium
KZ011-3082N060	6孔 6H	40	
KZ011-3082N070	7孔 7H	48	
KZ011-3082N080	8孔 8H	56	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws			

指骨锁定板 Phalange Locking Plates

指骨L型锁定板 Phalange L-type locking plate




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3072L050	5孔 5H	32	纯钛 Pure Titanium
KZ011-3072L060	6孔 6H	40	
KZ011-3072L070	7孔 7H	48	
KZ011-3072L080	8孔 8H	56	
KZ011-3072R050	5孔 5H	32	
KZ011-3072R060	6孔 6H	40	
KZ011-3072R070	7孔 7H	48	
KZ011-3072R080	8孔 8H	56	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws			

指骨斜L型锁定板 Phalange Slanted L-type locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3062L050	5孔 5H	34.5	纯钛 Pure Titanium
KZ011-3062L060	6孔 6H	42.5	
KZ011-3062L070	7孔 7H	50.5	
KZ011-3062L080	8孔 8H	58.5	
KZ011-3062R050	5孔 5H	34.5	
KZ011-3062R060	6孔 6H	42.5	
KZ011-3062R070	7孔 7H	50.5	
KZ011-3062R080	8孔 8H	58.5	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws			

1.5T型锁定接骨板 1.5 T-Type locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.071.200	4孔 4H	24.5	纯钛 Pure Titanium
101.071.201	5孔 5H	29.5	
101.071.202	6孔 6H	34.5	
101.071.203	7孔 7H	39.5	
101.071.204	8孔 8H	44.5	
备注：选用1.5锁定螺钉，1.5皮质骨螺钉 Remark:Select 1.5 locking screws,1.5 cortical screws			

1.5髌锁定接骨板 1.5 condylar locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.072.200	4孔 4H	26	纯钛 Pure Titanium
	101.072.201	5孔 5H	31	
	101.072.202	6孔 6H	36	
	备注：选用1.5锁定螺钉，1.5皮质骨螺钉 Remark:Select 1.5 locking screws,1.5 cortical screws			

指骨锁定板 Phalange Locking Plates

2.0直型锁定接骨板 2.0 Straight Locking Plate




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2012N030	3孔 3H	21	纯钛 Pure Titanium
KZ011-2012N040	4孔 4H	31	
KZ011-2012N050	5孔 5H	38	
KZ011-2012N060	6孔 6H	45	
KZ011-2012N070	7孔 7H	52	
KZ011-2012N080	8孔 8H	59	
备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨直型 II 型锁定板 Phalange Straight Type II Locking Plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1992N030	3孔 3H	20	纯钛 Pure Titanium
KZ011-1992N040	4孔 4H	27	
KZ011-1992N050	5孔 5H	34	
KZ011-1992N060	6孔 6H	41	
KZ011-1992N070	7孔 7H	48	
KZ011-1992N080	8孔 8H	55	
备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨Y型 II 型锁定板 Phalange Y-type II-type locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2002N030	3孔 3H	27	纯钛 Pure Titanium
KZ011-2002N040	4孔 4H	34	
KZ011-2002N050	5孔 5H	41	
KZ011-2002N060	6孔 6H	48	
KZ011-2002N070	7孔 7H	55	
备注：选用2.0锁定螺钉，2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			

指骨直型III型锁定板 Phalange Straight Type III Locking Plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2122N030	3孔 3H	21	纯钛 Pure Titanium
KZ011-2122N040	4孔 4H	29	
KZ011-2122N050	5孔 5H	37	
KZ011-2122N060	6孔 6H	45	
KZ011-2122N070	7孔 7H	53	
KZ011-2122N080	8孔 8H	61	
备注：选用螺钉HC2.0 Remark:Select HC2.0screws			


指骨锁定板 Phalange Locking Plates

指骨T型III型锁定板 Phalange T-type -III locking plates




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2032N040	4孔 4H	37	纯钛 Pure Titanium
KZ011-2032N050	5孔 5H	45	
KZ011-2032N060	6孔 6H	53	
KZ011-2032N070	7孔 7H	61	
KZ011-2032N090	8孔 8H	69	
备注：选用螺钉HC2.0 Remark:Select HC2.0screws			

指骨L型 I 型锁定板（左右式） L-Shaped Phalange I-type locking plates (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2042L040	4孔 4H	37	纯钛 Pure Titanium
KZ011-2042L050	5孔 5H	45	
KZ011-2042L060	6孔 6H	53	
KZ011-2042L070	7孔 7H	61	
KZ011-2042L080	8孔 8H	69	
KZ011-2042R040	4孔 4H	37	
KZ011-2042R050	5孔 5H	45	
KZ011-2042R060	6孔 6H	53	
KZ011-2042R070	7孔 7H	61	
KZ011-2042R080	8孔 8H	69	
备注：选用螺钉HC2.0 Remark:Select HC2.0 screws			

指骨L型锁定板120°（左右式） Curved 120° Phalange locking plates (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2052L040	4孔 4H	40.5	纯钛 Pure Titanium
KZ011-2052L050	5孔 5H	48.5	
KZ011-2052L060	6孔 6H	56.5	
KZ011-2052L070	7孔 7H	64.5	
KZ011-2052L080	8孔 8H	72.5	
KZ011-2052R040	4孔 4H	40.5	
KZ011-2052R050	5孔 5H	48.5	
KZ011-2052R060	6孔 6H	56.5	
KZ011-2052R070	7孔 7H	64.5	
KZ011-2052R080	8孔 8H	72.5	
备注：选用螺钉HC2.0 Remark:Select HC2.0 screws			

掌骨锁定板 Metacarpal locking plate

掌骨直型III型锁定板 Straight Metacarpals locking plate Type III

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2062N030	3孔 3H	25	纯钛 Pure Titanium
	KZ011-2062N040	4孔 4H	33	
	KZ011-2062N050	5孔 5H	41	
	KZ011-2062N060	6孔 6H	49	
	KZ011-2062N070	7孔 7H	57	
	KZ011-2062N080	7孔 8H	65	
备注: 选用2.7锁定螺钉, 2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws				

掌骨T型锁定板 I 型 T-Metacarpals locking plate Type I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2072N030	3孔 3H	32	纯钛 Pure Titanium
	KZ011-2072N040	4孔 4H	40	
	KZ011-2072N050	5孔 5H	48	
	KZ011-2072N060	6孔 6H	56	
备注: 选用2.7锁定螺钉, 2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws				

掌骨L型锁定板 I 型（左右式） T-Metacarpals locking plate Type I (left and right type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2082L030	3孔 3H	32	纯钛 Pure Titanium
	KZ011-2082L040	4孔 4H	40	
	KZ011-2082L050	5孔 5H	48	
	KZ011-2082L060	6孔 6H	56	
	KZ011-2082R030	3孔 3H	32	
	KZ011-2082R040	4孔 4H	40	
	KZ011-2082R050	5孔 5H	48	
	KZ011-2082R060	6孔 6H	56	
备注: 选用2.7锁定螺钉, 2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws				

掌骨T型锁定板 III 型 T-Metacarpals locking plate Type III

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2132N040	4孔 4H	39	纯钛 Pure Titanium
	KZ011-2132N050	5孔 5H	47	
	KZ011-2132N060	6孔 6H	55	
	KZ011-2132N070	7孔 7H	63	
	KZ011-2132N080	8孔 8H	71	
备注: 选用螺钉HC2.7 Remark:Select HC2.7screws				

掌骨锁定板 Metacarpal locking plate

掌骨L型锁定板 II 型（左右式） L-Metacarpals locking plate Type II (left and right type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2092L030	3孔 3H	35	纯钛 Pure Titanium
	KZ011-2092L040	4孔 4H	43	
	KZ011-2092L050	5孔 5H	51	
	KZ011-2092L060	6孔 6H	59	
	KZ011-2092R030	3孔 3H	35	
	KZ011-2092R040	4孔 4H	43	
	KZ011-2092R050	5孔 5H	51	
	KZ011-2092R060	6孔 6H	59	
备注: 选用2.7锁定螺钉, 2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws				

掌骨Y型锁定板 I 型 Y-Metacarpals locking plate Type I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2102N030	3孔 3H	34	纯钛 Pure Titanium
	KZ011-2102N040	4孔 4H	42	
	KZ011-2102N050	5孔 5H	50	
	KZ011-2102N060	6孔 6H	58	
	KZ011-2102N070	7孔 7H	66	
备注: 选用2.7锁定螺钉, 2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws				

掌骨直型 I 型锁定板 Straight Metacarpals locking plate Type I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2112N030	3孔 3H	24	纯钛 Pure Titanium
	KZ011-2112N040	4孔 4H	32	
	KZ011-2112N050	5孔 5H	40	
	KZ011-2112N060	6孔 6H	48	
	KZ011-2112N070	7孔 7H	56	
	KZ011-2112N080	8孔 8H	64	
备注: 选用螺钉HC2.7 Remark:Select HC2.7 screws				

掌骨锁定板Metacarpal locking plate

2.4T型锁定接骨板2.4 T type locking bone plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.080.201	4孔 4H	36.25	纯钛 Pure Titanium
	101.080.202	5孔 5H	43.5	
	101.080.203	6孔 6H	50.75	
	101.080.204	7孔 7H	58	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4髁锁定接骨板2.4 Condylar locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.081.201	4孔 4H	37.25	纯钛 Pure Titanium
	101.081.202	5孔 5H	44.5	
	101.081.203	6孔 6H	51.75	
	101.081.204	7孔 7H	59	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4Y型锁定接骨板2.4 Y type locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.082.201	4孔 4H	38.25	纯钛 Pure Titanium
	101.082.202	5孔 5H	45.5	
	101.082.203	6孔 6H	52.75	
	101.082.204	7孔 7H	60	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4T型锁定接骨板II2.4 T type locking bone plate II

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.083.201	4孔 4H	36.25	纯钛 Pure Titanium
	101.083.202	5孔 5H	43.5	
	101.083.203	6孔 6H	50.75	
	101.083.204	7孔 7H	58	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

掌骨锁定板Metacarpal locking plate

2.7髁锁定接骨板2.7 Condylar locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.087.201	4孔 4H	41	纯钛 Pure Titanium
	101.087.202	5孔 5H	49	
	101.087.203	6孔 6H	57	
	101.087.204	7孔 7H	65	
	备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

掌骨L型I型锁定板（左右式）L-Metacarpals locking plate Type I (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2142L040	4孔 4H	39	纯钛 Pure Titanium
KZ011-2142L050	5孔 5H	47	
KZ011-2142L060	6孔 6H	55	
KZ011-2142L070	7孔 7H	63	
KZ011-2142L080	8孔 8H	71	
KZ011-2142R040	4孔 4H	39	
KZ011-2142R050	5孔 5H	47	
KZ011-2142R060	6孔 6H	55	
KZ011-2142R070	7孔 7H	63	
KZ011-2142R080	8孔 8H	71	
备注：选用螺钉HC2.7 Remark:Select HC2.7screws			

掌骨L型锁定板120°（左右式）Curved 120° Metacarpals locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2152L040	4孔 4H	42.5	纯钛 Pure Titanium
KZ011-2152L050	5孔 5H	50.5	
KZ011-2152L060	6孔 6H	58.5	
KZ011-2152L070	7孔 7H	66.5	
KZ011-2152L080	8孔 8H	74.5	
KZ011-2152R040	4孔 4H	42.5	
KZ011-2152R050	5孔 5H	50.5	
KZ011-2152R060	6孔 6H	58.5	
KZ011-2152R070	7孔 7H	66.5	
KZ011-2152R080	8孔 8H	74.5	
备注：选用螺钉HC2.7 Remark:Select HC2.7screws			

上肢锁定板Upper limb locking plate

锁骨钩锁定板- IIClavicle hook locking plate- II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1322L041	4孔 4H	70	纯钛 Pure Titanium
KZ011-1322L051	5孔 5H	85	
KZ011-1322L061	6孔 6H	99	
KZ011-1322L071	7孔 7H	112	
KZ011-1322R041	4孔 4H	70	
KZ011-1322R051	5孔 5H	85	
KZ011-1322R061	6孔 6H	99	
KZ011-1322R071	7孔 7H	112	
KZ011-1322L042	4孔 4H	70	
KZ011-1322L052	5孔 5H	85	
KZ011-1322L062	6孔 6H	99	
KZ011-1322L072	7孔 7H	112	
KZ011-1322R042	4孔 4H	70	
KZ011-1322R052	5孔 5H	85	
KZ011-1322R062	6孔 6H	99	
KZ011-1322R072	7孔 7H	112	
KZ011-1322L043	4孔 4H	70	
KZ011-1322L053	5孔 5H	85	
KZ011-1322L063	6孔 6H	99	
KZ011-1322L073	7孔 7H	112	
KZ011-1322R043	4孔 4H	70	
KZ011-1322R053	5孔 5H	85	
KZ011-1322R063	6孔 6H	99	
KZ011-1322R073	7孔 7H	112	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

3.5锁骨锁定板 I3.5Clavicle locking plate I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1022N030	3孔 3H	94	纯钛 Pure Titanium
KZ011-1022N040	4孔 4H	107	
KZ011-1022N050	5孔 5H	120	
KZ011-1022N060	6孔 6H	133	
KZ011-1022N070	7孔 7H	146	
KZ011-1022N080	8孔 8H	159	
KZ011-1022N090	9孔 9H	172	
KZ011-1022N100	10孔 10H	185	
KZ011-1022N110	11孔 11H	198	
KZ011-1022N120	12孔 12H	211	
KZ011-1022N130	13孔 13H	224	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

上肢锁定板Upper limb locking plate

锁骨锁定板Clavicle locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1452L031	3孔 3H	43	纯钛 Pure Titanium
KZ011-1452L041	4孔 4H	57.75	
KZ011-1452L051	5孔 5H	72.5	
KZ011-1452L061	6孔 6H	87.25	
KZ011-1452L071	7孔 7H	102	
KZ011-1452L081	8孔 8H	116.75	
KZ011-1452L091	9孔 9H	131.5	
KZ011-1452L101	10孔 10H	146.25	
KZ011-1452L111	11孔 11H	161	
KZ011-1452L121	12孔 12H	175.75	
KZ011-1452R031	3孔 3H	43	
KZ011-1452R041	4孔 4H	57.75	
KZ011-1452R051	5孔 5H	72.5	
KZ011-1452R061	6孔 6H	87.25	
KZ011-1452R071	7孔 7H	102	
KZ011-1452R081	8孔 8H	116.75	
KZ011-1452R091	9孔 9H	131.5	
KZ011-1452R101	10孔 10H	146.25	
KZ011-1452R111	11孔 11H	161	
KZ011-1452R121	12孔 12H	175.75	
备注：选用3.5锁定螺钉，4.0皮质骨螺钉 Remark:Select 3.5 locking screws,4.0 cortical screws			

锁骨前侧型锁定板Anterior clavicle locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1412L030	3孔 3H	69	纯钛 Pure Titanium
KZ011-1412L040	4孔 4H	81	
KZ011-1412L050	5孔 5H	94	
KZ011-1412L060	6孔 6H	108	
KZ011-1412L070	7孔 7H	123	
KZ011-1412L080	8孔 8H	135	
KZ011-1412L090	9孔 9H	145	
KZ011-1412R030	3孔 3H	69	
KZ011-1412R040	4孔 4H	81	
KZ011-1412R050	5孔 5H	94	
KZ011-1412R060	6孔 6H	108	
KZ011-1412R070	7孔 7H	123	
KZ011-1412R080	8孔 8H	135	
KZ011-1412R090	9孔 9H	145	
备注：选用2.4/2.7锁定接骨螺钉，3.5锁定螺钉，4.0皮质骨螺钉 Remark:Select 2.4/2.7 locking bone screws, 3.5 locking screws, and 4.0 cortical screws.			

上肢锁定板Upper limb locking plate

锁骨内侧上方锁定板Clavicle Inner upper locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1392L040	4孔 4H	46.4	纯钛 Pure Titanium
KZ011-1392L050	5孔 5H	57.5	
KZ011-1392L060	6孔 6H	68.4	
KZ011-1392L070	7孔 7H	80	
KZ011-1392L080	8孔 8H	97	
KZ011-1392L090	9孔 9H	109	
KZ011-1392L100	10孔 10H	121	
KZ011-1392L110	11孔 11H	133	
KZ011-1392L120	12孔 12H	145	
KZ011-1392R040	4孔 4H	46.4	
KZ011-1392R050	5孔 5H	57.5	
KZ011-1392R060	6孔 6H	68.4	
KZ011-1392R070	7孔 7H	80	
KZ011-1392R080	8孔 8H	97	
KZ011-1392R090	9孔 9H	109	
KZ011-1392R100	10孔 10H	121	
KZ011-1392R110	11孔 11H	133	
KZ011-1392R120	12孔 12H	145	
备注：选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

锁骨干重建型锁定板Clavicle trunk reconstruction locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1082L050	5孔 5H	74	纯钛 Pure Titanium
KZ011-1082L060	6孔 6H	99	
KZ011-1082L070	7孔 7H	109	
KZ011-1082L080	8孔 8H	121	
KZ011-1082L090	9孔 9H	124.8	
KZ011-1082L100	10孔 10H	136.4	
KZ011-1082R050	5孔 5H	74	
KZ011-1082R060	6孔 6H	99	
KZ011-1082R070	7孔 7H	109	
KZ011-1082R080	8孔 8H	121	
KZ011-1082R090	9孔 9H	124.8	
KZ011-1082R100	10孔 10H	136.4	
备注：选用3.5锁定螺钉， 3.5皮质骨螺钉,4.0松质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws, and 4.0 cancellous screws.			

上肢锁定板Upper limb locking plate

2.7直型锁定板2.7 Straight Locking Plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.098.200	4孔 4H	40	纯钛 Pure Titanium
101.098.201	5孔 5H	49	
101.098.202	6孔 6H	58	
101.098.203	7孔 7H	67	
101.098.204	8孔 8H	76	
101.098.205	9孔 9H	85	
101.098.206	10孔 10H	97	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

2.0重建锁定板2.0 Rebuild Locking Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.103.200	6孔 6H	40.5	纯钛 Pure Titanium
	101.103.201	12孔 12H	81	
	备注: 选用2.0锁定螺钉, 2.0皮质骨螺钉 Remark: Select 2.0 locking screws, 2.0 cortical screws			

2.4重建锁定板2.4 Rebuild Locking Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.119.200	6孔 6H	44	纯钛 Pure Titanium
	101.119.201	12孔 12H	88	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

3.5重建锁定板 I3.5 Rebuild Locking Plate I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.105.200	3孔 3H	39	纯钛 Pure Titanium
101.105.202	5孔 5H	65	
101.105.204	7孔 7H	91	
101.105.206	9孔 9H	117	
101.105.208	11孔 11H	143	
101.105.210	13孔 13H	169	
101.105.212	15孔 15H	195	
101.105.215	18孔 18H	234	
备注：选用螺钉HC3.5 Remark:Select screws HC3.5			

上肢锁定板Upper limb locking plate

管型锁定板 I Tubular locking plate I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1952N030	3孔 3H	38	纯钛 Pure Titanium
KZ011-1952N040	4孔 4H	50	
KZ011-1952N050	5孔 5H	62	
KZ011-1952N060	6孔 6H	74	
KZ011-1952N070	7孔 7H	86	
KZ011-1952N080	8孔 8H	98	
KZ011-1952N090	9孔 9H	110	
KZ011-1952N100	10孔 10H	122	
KZ011-1952N110	11孔 11H	134	
KZ011-1952N120	12孔 12H	146	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

上肢不等厚干骺端锁定板Upper limb unequal thickness epiphyseal locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1202N050	5孔 5H	73	纯钛 Pure Titanium
KZ011-1202N060	6孔 6H	86	
KZ011-1202N070	7孔 7H	99	
KZ011-1202N080	8孔 8H	112	
KZ011-1202N090	9孔 9H	125	
KZ011-1202N100	10孔 10H	138	
KZ011-1202N110	11孔 11H	151	
KZ011-1202N120	12孔 12H	164	
KZ011-1202N140	14孔 14H	190	
KZ011-1202N160	16孔 16H	216	
KZ011-1202N180	18孔 18H	242	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

有限接触锁定板Limited contact locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1382N040	4孔 4H	59	纯钛 Pure Titanium
KZ011-1382N050	5孔 5H	72	
KZ011-1382N060	6孔 6H	85	
KZ011-1382N070	7孔 7H	98	
KZ011-1382N080	8孔 8H	111	
KZ011-1382N090	9孔 9H	124	
KZ011-1382N100	10孔 10H	137	
KZ011-1382N110	11孔 11H	150	
KZ011-1382N120	12孔 12H	163	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉,4.0松质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws, and 4.0 cancellous screws.			

上肢锁定板Upper limb locking plate

肱骨远端内侧锁定板Distal medial locking plate of humerus



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1362L030	3孔 3H	58	纯钛 Pure Titanium
KZ011-1362L040	4孔 4H	71	
KZ011-1362L050	5孔 5H	84	
KZ011-1362L060	6孔 6H	97	
KZ011-1362L070	7孔 7H	110	
KZ011-1362L080	8孔 8H	123	
KZ011-1362L090	9孔 9H	136	
KZ011-1362L100	10孔 10H	149	
KZ011-1362L110	11孔 11H	162	
KZ011-1362L120	12孔 12H	175	
KZ011-1362L130	13孔 13H	188	
KZ011-1362L140	14孔 14H	201	
KZ011-1362R030	3孔 3H	58	
KZ011-1362R040	4孔 4H	71	
KZ011-1362R050	5孔 5H	84	
KZ011-1362R060	6孔 6H	97	
KZ011-1362R070	7孔 7H	110	
KZ011-1362R080	8孔 8H	123	
KZ011-1362R090	9孔 9H	136	
KZ011-1362R100	10孔 10H	149	
KZ011-1362R110	11孔 11H	162	
KZ011-1362R120	12孔 12H	175	
KZ011-1362R130	13孔 13H	188	
KZ011-1362R140	14孔 14H	201	
备注：选用2.4/2.7/3.5锁定螺钉，4.0皮质骨螺钉 Remark:Select 2.4/2.7/3.5 locking screws and 4.0 cortical screws.			

肱骨远端背外侧带支撑锁定板Distal dorsolateral humeral band support locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1352L030	3孔 3H	65	纯钛 Pure Titanium
KZ011-1352L040	4孔 4H	78	
KZ011-1352L050	5孔 5H	91	
KZ011-1352L060	6孔 6H	104	
KZ011-1352L070	7孔 7H	117	
KZ011-1352L080	8孔 8H	130	
KZ011-1352L090	9孔 9H	143	
KZ011-1352L110	11孔 11H	169	
KZ011-1352L140	14孔 14H	208	
KZ011-1352R030	3孔 3H	65	
KZ011-1352R040	4孔 4H	78	
KZ011-1352R050	5孔 5H	91	
KZ011-1352R060	6孔 6H	104	
KZ011-1352R070	7孔 7H	117	
KZ011-1352R080	8孔 8H	130	
KZ011-1352R090	9孔 9H	143	
KZ011-1352R110	11孔 11H	169	
KZ011-1352R140	14孔 14H	208	
备注：选用2.4/2.7/3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 2.4/2.7/3.5 locking screws and 3.5 cortical screws.			

上肢锁定板Upper limb locking plate

尺骨鹰嘴锁定板Ulnar olecranon locking plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1112L020	2孔 2H	86	纯钛 Pure Titanium
KZ011-1112L030	3孔 3H	103	
KZ011-1112L040	4孔 4H	112	
KZ011-1112L050	5孔 5H	129	
KZ011-1112L060	6孔 6H	138	
KZ011-1112L070	7孔 7H	155	
KZ011-1112L080	8孔 8H	164	
KZ011-1112L090	9孔 9H	181	
KZ011-1112L100	10孔 10H	190	
KZ011-1112L110	11孔 11H	201	
KZ011-1112L120	12孔 12H	216	
KZ011-1112R020	2孔 2H	86	
KZ011-1112R030	3孔 3H	103	
KZ011-1112R040	4孔 4H	112	
KZ011-1112R050	5孔 5H	129	
KZ011-1112R060	6孔 6H	138	
KZ011-1112R070	7孔 7H	155	
KZ011-1112R080	8孔 8H	164	
KZ011-1112R090	9孔 9H	181	
KZ011-1112R100	10孔 10H	190	
KZ011-1112R110	11孔 11H	201	
KZ011-1112R120	12孔 12H	216	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉,4.0松质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws, and 4.0 cancellous screws.			

尺桡骨直型锁定板-1Straight locking plate for radius and ulna-1



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1282N040	4孔 4H	59	纯钛 Pure Titanium
KZ011-1282N050	5孔 5H	72	
KZ011-1282N060	6孔 6H	85	
KZ011-1282N070	7孔 7H	98	
KZ011-1282N080	8孔 8H	111	
KZ011-1282N090	9孔 9H	124	
KZ011-1282N100	10孔 10H	137	
KZ011-1282N110	11孔 11H	150	
KZ011-1282N120	12孔 12H	163	
备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark: Select 3.5 locking screws, 3.5 cortical screws			

重建型直型锁定板Reconstruction type straight locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1262N030	3孔 3H	36	
KZ011-1262N040	4孔 4H	48	
KZ011-1262N050	5孔 5H	58	
KZ011-1262N060	6孔 6H	70	
KZ011-1262N070	7孔 7H	82	
KZ011-1262N080	8孔 8H	94	
KZ011-1262N090	9孔 9H	106	
KZ011-1262N100	10孔 10H	118	

上肢锁定板Upper limb locking plate

				纯钛 Pure Titanium
	KZ011-1262N110	11孔 11H	130	
	KZ011-1262N120	12孔 12H	142	
	KZ011-1262N130	13孔 13H	154	
	KZ011-1262N140	14孔 14H	166	
	KZ011-1262N150	15孔 15H	178	
	KZ011-1262N160	16孔 16H	190	
	KZ011-1262N180	18孔 18H	214	
	KZ011-1262N200	20孔 20H	238	
	KZ011-1262N220	22孔 22H	262	
	备注：选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

重建弧型锁定板Reconstruct the arc-shaped locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1212N040	4孔 4H	46	纯钛 Pure Titanium
	KZ011-1212N050	5孔 5H	58	
	KZ011-1212N060	6孔 6H	70	
	KZ011-1212N070	7孔 7H	82	
	KZ011-1212N080	8孔 8H	94	
	KZ011-1212N090	9孔 9H	106	
	KZ011-1212N100	10孔 10H	118	
	KZ011-1212N110	11孔 11H	130	
	KZ011-1212N120	12孔 12H	142	
	KZ011-1212N130	13孔 13H	154	
	KZ011-1212N140	14孔 14H	166	
	备注：选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

尺桡骨锁定板（圆孔型）Radius and ulna locking plate (round hole type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2262N050	5孔 5H	72.2	纯钛 Pure Titanium
	KZ011-2262N060	6孔 6H	85.2	
	KZ011-2262N070	7孔 7H	98.2	
	KZ011-2262N080	8孔 8H	111.2	
	KZ011-2262N090	9孔 9H	124.2	
	KZ011-2262N100	10孔 10H	137.2	
	备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

重建型儿童锁定板Reconstructed Child Locking Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2272N040	4孔 4H	40	纯钛 Pure Titanium
	KZ011-2272N050	5孔 5H	50	
	KZ011-2272N060	6孔 6H	60	
	KZ011-2272N070	7孔 7H	70	
	KZ011-2272N080	8孔 8H	80	
	KZ011-2272N100	10孔 10H	90	
	备注：选用螺钉HC2.7 Remark:Select screws HC2.7			

上肢锁定板Upper limb locking plate

1/3圆锁定板1/3 circle locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2282N040	4孔 4H	49	纯钛 Pure Titanium
KZ011-2282N050	5孔 5H	61	
KZ011-2282N060	6孔 6H	73	
KZ011-2282N070	7孔 7H	85	
KZ011-2282N080	8孔 8H	97	
KZ011-2282N090	9孔 9H	109	
KZ011-2282N100	10孔 10H	121	
备注：选用螺钉HC3.5 Remark:Select screws HC3.5			

鹰嘴锁定板 II Eagle Beak Locking Plate II

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2292N080	8孔 8H	46	纯钛 Pure Titanium
	备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

尺骨鹰嘴锁定板IIIUlnar olecranon locking plate III



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2302N030	3孔 3H	55	纯钛 Pure Titanium
KZ011-2302N040	4孔 4H	69	
KZ011-2302N050	5孔 5H	83	
KZ011-2302N060	6孔 6H	97	
KZ011-2302N070	7孔 7H	111	
KZ011-2302N080	8孔 8H	125	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

桡骨锁定板Radial locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2312N030	3孔 3H	70.5	纯钛 Pure Titanium
KZ011-2312N040	4孔 4H	82.5	
KZ011-2312N050	5孔 5H	94.5	
KZ011-2312N060	6孔 6H	106.5	
KZ011-2312N070	7孔 7H	118.5	
备注：选用螺钉HC3.5 Remark:Select screws HC3.5			

上肢锁定板Upper limb locking plate

肱骨近端 II 型锁定板（左右式）Proximal humerus Type II locking plate (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2322L050	/	43.3	纯钛 Pure Titanium
	KZ011-2322R050	/	43.3	
	备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

肱骨近端髁锁定板IV型Proximal humeral condyle locking plate Type IV



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2332N040	4孔 4H	74	纯钛 Pure Titanium
KZ011-2332N050	5孔 5H	86	
KZ011-2332N060	6孔 6H	98	
KZ011-2332N070	7孔 7H	110	
KZ011-2332N080	8孔 8H	122	
KZ011-2332N090	9孔 9H	134	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

肱骨颈（左右式）Humeral neck (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2342L030	3孔 3H	85	纯钛 Pure Titanium
KZ011-2342L040	4孔 4H	98	
KZ011-2342L050	5孔 5H	111	
KZ011-2342L060	6孔 6H	124	
KZ011-2342L070	7孔 7H	137	
KZ011-2342L080	8孔 8H	150	
KZ011-2342L090	9孔 9H	163	
KZ011-2342R030	3孔 3H	85	
KZ011-2342R040	4孔 4H	98	
KZ011-2342R050	5孔 5H	111	
KZ011-2342R060	6孔 6H	124	
KZ011-2342R070	7孔 7H	137	
KZ011-2342R080	8孔 8H	150	
KZ011-2342R090	9孔 9H	163	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

肱骨远端内侧 II 型Distal medial humerus type II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2362N060	6孔 6H	112	纯钛 Pure Titanium
KZ011-2362N080	8孔 8H	138	
KZ011-2362N100	10孔 10H	164	
KZ011-2362N120	12孔 12H	190	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

上肢锁定板Upper limb locking plate

肱骨远端外侧锁定板（左右式）Lateral locking plate of distal humerus (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2352L040	4孔 4H	96.6	纯钛 Pure Titanium
KZ011-2352L060	6孔 6H	129.6	
KZ011-2352L070	7孔 7H	143.6	
KZ011-2352L080	8孔 8H	161.6	
KZ011-2352L100	10孔 10H	189.6	
KZ011-2352R040	4孔 4H	96.6	
KZ011-2352R060	6孔 6H	129.6	
KZ011-2352R070	7孔 7H	143.6	
KZ011-2352R080	8孔 8H	161.6	
KZ011-2352R100	10孔 10H	189.6	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

肱骨远端Y型锁定板（左右式）Y-shaped locking plate for distal humerus (left and right).



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2372L030	3孔 3H	70	纯钛 Pure Titanium
KZ011-2372L040	4孔 4H	82	
KZ011-2372L050	5孔 5H	94	
KZ011-2372L060	6孔 6H	106	
KZ011-2372L070	7孔 7H	118	
KZ011-2372L080	8孔 8H	130	
KZ011-2372R030	3孔 3H	70	
KZ011-2372R040	4孔 4H	82	
KZ011-2372R050	5孔 5H	94	
KZ011-2372R060	6孔 6H	106	
KZ011-2372R070	7孔 7H	118	
KZ011-2372R080	8孔 8H	130	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

S型锁骨 II 型锁定板（左右式）S-shaped clavicle type II locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2402L040	4孔 4H	46	纯钛 Pure Titanium
KZ011-2402L050	5孔 5H	58	
KZ011-2402L060	6孔 6H	70	
KZ011-2402L070	7孔 7H	82	
KZ011-2402L080	8孔 8H	93	
KZ011-2402L090	9孔 9H	104	
KZ011-2402L100	10孔 10H	115	
KZ011-2402L120	12孔 12H	135	
KZ011-2402R040	4孔 4H	46	
KZ011-2402R050	5孔 5H	58	
KZ011-2402R060	6孔 6H	70	
KZ011-2402R070	7孔 7H	82	
KZ011-2402R080	8孔 8H	93	
KZ011-2402R090	9孔 9H	104	
KZ011-2402R100	10孔 10H	115	
KZ011-2402R120	12孔 12H	135	
备注：选用螺钉HC3.5 Remark:Select screws HC3.5			

上肢锁定板Upper limb locking plate

锁骨远端I型（左右式）组合孔Distal clavicle type I (left and right) combination hole



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2382L040	4孔 4H	49.3	纯钛 Pure Titanium
KZ011-2382L050	5孔 5H	60.2	
KZ011-2382L060	6孔 6H	71.3	
KZ011-2382L070	7孔 7H	82.9	
KZ011-2382L080	8孔 8H	94.8	
KZ011-2382L090	9孔 9H	106.8	
KZ011-2382L100	10孔 10H	118.7	
KZ011-2382L110	11孔 11H	130.1	
KZ011-2382L120	12孔 12H	141	
KZ011-2382R040	4孔 4H	49.3	
KZ011-2382R050	5孔 5H	60.2	
KZ011-2382R060	6孔 6H	71.3	
KZ011-2382R070	7孔 7H	82.9	
KZ011-2382R080	8孔 8H	94.8	
KZ011-2382R090	9孔 9H	106.8	
KZ011-2382R100	10孔 10H	118.7	
KZ011-2382R110	11孔 11H	130.1	
KZ011-2382R120	12孔 12H	141	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

锁骨远端I型（左右式）Distal clavicle type I (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2392L040	4孔 4H	46.4	纯钛 Pure Titanium
KZ011-2392L050	5孔 5H	57.5	
KZ011-2392L060	6孔 6H	68.4	
KZ011-2392L070	7孔 7H	80	
KZ011-2392L080	8孔 8H	91.9	
KZ011-2392L090	9孔 9H	103.9	
KZ011-2392L100	10孔 10H	115.7	
KZ011-2392L110	11孔 11H	127.2	
KZ011-2392L120	12孔 12H	138	
KZ011-2392R040	4孔 4H	46.4	
KZ011-2392R050	5孔 5H	57.5	
KZ011-2392R060	6孔 6H	68.4	
KZ011-2392R070	7孔 7H	80	
KZ011-2392R080	8孔 8H	91.9	
KZ011-2392R090	9孔 9H	103.9	
KZ011-2392R100	10孔 10H	115.7	
KZ011-2392R110	11孔 11H	127.2	
KZ011-2392R120	12孔 12H	138	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

上肢锁定板 Upper limb locking plate

S型锁骨 I 型锁定板（左右式） S-Type clavicle type I (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2432L040	4孔 4H	46	纯钛 Pure Titanium
KZ011-2432L050	5孔 5H	58	
KZ011-2432L060	6孔 6H	70	
KZ011-2432L070	7孔 7H	82	
KZ011-2432L080	8孔 8H	93	
KZ011-2432L090	9孔 9H	104	
KZ011-2432L100	10孔 10H	115	
KZ011-2432L120	12孔 12H	135	
KZ011-2432R040	4孔 4H	46	
KZ011-2432R050	5孔 5H	58	
KZ011-2432R060	6孔 6H	70	
KZ011-2432R070	7孔 7H	82	
KZ011-2432R080	8孔 8H	93	
KZ011-2432R090	9孔 9H	104	
KZ011-2432R100	10孔 10H	115	
KZ011-2432R120	12孔 12H	135	
备注：选用螺钉HC3.5 Remark:Select screws HC3.5			

三叶型锁定板（左右式） Trilobal locking plate (left and right type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2442L090	/	88.7	纯钛 Pure Titanium
	KZ011-2442R090	/	88.7	
	备注：选用螺钉HC2.7 Remark:Select screws HC2.7			

弓 I 型锁定板（左右式） Bow-type locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2452L080	/	77.4	纯钛 Pure Titanium
KZ011-2452R080	/	77.4	
备注：选用螺钉HC2.7 Remark:Select screws HC2.7			

成角锁定板（左右式） Angle locking plate (left and right type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2462L040	/	21	纯钛 Pure Titanium
	KZ011-2462R040	/	21	
	备注：选用螺钉HC2.7 Remark:Select screws HC2.7			

上肢锁定板 Upper limb locking plate

弓 II 型锁定板（左右式） Bow Type II Locking Plate (Left and Right Type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2472L100	/	84.9	纯钛 Pure Titanium
	KZ011-2472R100	/	84.9	
	备注：选用螺钉HC2.7 Remark:Select screws HC2.7			

重建锁定板（组合孔） Reconstruction locking plate (combination hole)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2492N030	3孔 3H	40.2	纯钛 Pure Titanium
KZ011-2492N040	4孔 4H	53.2	
KZ011-2492N050	5孔 5H	66.2	
KZ011-2492N060	6孔 6H	79.2	
KZ011-2492N070	7孔 7H	92.2	
KZ011-2492N080	8孔 8H	105.2	
KZ011-2492N090	9孔 9H	118.2	
KZ011-2492N100	10孔 10H	131.2	
KZ011-2492N110	11孔 11H	144.2	
KZ011-2492N120	12孔 12H	157.2	
KZ011-2492N130	13孔 13H	170.2	
KZ011-2492N140	14孔 14H	183.2	
KZ011-2492N150	15孔 15H	196.2	
KZ011-2492N160	16孔 16H	209.2	
备注：选用螺钉HC3.5.HC4.0 Remark:Select screws HC3.5.HC4.0			

骨盆锁定板 Pelvic locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2482N040	4孔 4H	48	纯钛 Pure Titanium
KZ011-2482N050	5孔 5H	56.45	
KZ011-2482N060	6孔 6H	64.9	
KZ011-2482N070	7孔 7H	73.35	
KZ011-2482N080	8孔 8H	81.8	
KZ011-2482N090	9孔 9H	90.25	
KZ011-2482N100	10孔 10H	98.7	
KZ011-2482N110	11孔 11H	107.15	
KZ011-2482N120	12孔 12H	115.6	
KZ011-2482N130	13孔 13H	124.05	
KZ011-2482N140	14孔 14H	132.5	
备注：选用螺钉HC3.5.HC4.0 Remark:Select screws HC3.5.HC4.0			


上肢锁定板Metacarpal locking plate

3.5锁骨锁定板III3.5 Clavicle locking plate III



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.005.200	3孔 3H	69	纯钛 Pure Titanium
101.005.202	4孔 4H	81	
101.005.204	5孔 5H	94	
101.005.206	6孔 6H	108	
101.005.208	7孔 7H	123	
101.005.210	8孔 8H	135	
101.005.201	3孔 3H	69	
101.005.203	4孔 4H	81	
101.005.205	5孔 5H	94	
101.005.207	6孔 6H	108	
101.005.209	7孔 7H	123	
101.005.211	8孔 8H	135	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

3.5肱骨远端外侧锁定板I3.5 Lateral locking plate of distal humerus I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.010.204	5孔 5H	91	纯 钛 Pure Titanium
101.010.208	7孔 7H	117	
101.010.212	9孔 9H	143	
101.010.216	11孔 11H	169	
101.010.205	5孔 5H	91	
101.010.209	7孔 7H	117	
101.010.213	9孔 9H	143	
101.010.217	11孔 11H	169	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

2.0尺骨茎突锁定板2.0 Ulnar styloid process locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.016.200	7孔 7H	46	纯 钛 Pure Titanium
备注：选用2.0锁定螺钉,2.0皮质骨螺钉 Remark:Select 2.0 locking screws,2.0 cortical screws			


2.4桡骨头锁定板II2.4 Radial head locking plate II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.031.200	2孔 2H	37.5	纯 钛 Pure Titanium
101.031.202	3孔 3H	46.5	
101.031.204	4孔 4H	55	
101.031.206	5孔 5H	64.5	
101.031.201	2孔 2H	37.5	
101.031.203	3孔 3H	46.5	
101.031.205	4孔 4H	55	
101.031.207	5孔 5H	64.5	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

上肢锁定板Metacarpal locking plate

3.5 J型重建锁定接骨板I3.5 J type reconstruction locking plate I




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.035.200	10孔 10H	113	纯 钛 Pure Titanium
101.035.202	10孔 11H	126	
101.035.204	10孔 12H	139	
101.035.208	10孔 14H	165	
101.035.212	10孔 16H	191	
101.035.201	10孔 10H	113	
101.035.203	10孔 11H	126	
101.035.205	10孔 12H	139	
101.035.209	10孔 14H	165	
101.035.213	10孔 16H	191	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

肱骨近端锁定板-5Proximal humeral locking plate-5



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1022N030	3孔 3H	94	纯钛 Pure Titanium
KZ011-1022N040	4孔 4H	107	
KZ011-1022N050	5孔 5H	120	
KZ011-1022N060	6孔 6H	133	
KZ011-1022N070	7孔 7H	146	
KZ011-1022N080	8孔 8H	159	
KZ011-1022N090	9孔 9H	172	
KZ011-1022N100	10孔 10H	185	
KZ011-1022N110	11孔 11H	198	
KZ011-1022N120	12孔 12H	211	
KZ011-1022N130	13孔 13H	224	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

2.4直型锁定板2.4 Straight Locking Plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.097.200	4孔 4H	36	纯钛 Pure Titanium
101.097.201	5孔 5H	44	
101.097.202	6孔 6H	52	
101.097.203	7孔 7H	60	
101.097.204	8孔 8H	68	
备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

桡骨远端锁定板

Distal radius locking plate

2.4斜L型桡骨远端背侧锁定板 I

2.4 Oblique L-shaped dorsal locking plate of distal radius I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.020.200	3孔 3H	43	纯钛 Pure Titanium
	101.020.202	4孔 4H	52	
	101.020.204	5孔 5H	61	
	101.020.201	3孔 3H	43	
	101.020.203	4孔 4H	52	
	101.020.205	5孔 5H	61	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.				

2.4斜L型桡骨远端背侧锁定板 II

2.4 Oblique L-shaped dorsal locking plate II for distal radius

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.021.200	3孔 3H	42	纯钛 Pure Titanium
	101.021.202	4孔 4H	51	
	101.021.204	5孔 5H	60	
	101.021.201	3孔 3H	42	
	101.021.203	4孔 4H	51	
	101.021.205	5孔 5H	60	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.				

2.4L型桡骨远端背侧锁定板 II

2.4 L type distal radius dorsal locking plate II

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.019.200	3孔 3H	40	纯钛 Pure Titanium
	101.019.202	4孔 4H	49	
	101.019.204	5孔 5H	58	
	101.019.201	3孔 3H	40	
	101.019.203	4孔 4H	49	
	101.019.205	5孔 5H	58	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.				

2.4L型桡骨远端背侧锁定板 I

2.4 L type distal radius dorsal locking plate I

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.018.200	3孔 3H	40	纯钛 Pure Titanium
	101.018.202	4孔 4H	49	
	101.018.204	5孔 5H	58	
	101.018.201	3孔 3H	40	
	101.018.203	4孔 4H	49	
	101.018.205	5孔 5H	58	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.				

桡骨远端锁定板

Distal radius locking plate


桡骨远端T型锁定板

Distal radius T-locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1532N030	3孔 3H	40	纯钛 Pure Titanium
	KZ011-1532N040	4孔 4H	49	
	KZ011-1532N050	5孔 5H	58	
	KZ011-1532N060	6孔 6H	67	
	KZ011-1532N070	7孔 7H	76	
	备注：选用2.7锁定螺钉,2.7皮质骨螺钉 Remark:Select 2.7 locking screws,2.7 cortical screws			

桡骨远端T型弯锁定板

Distal radius T-curve locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1562N031	3孔 3H	45	纯钛 Pure Titanium
	KZ011-1562N041	4孔 4H	56	
	KZ011-1562N051	5孔 5H	67	
	KZ011-1562N061	6孔 6H	78	
	KZ011-1562N071	7孔 7H	89	
	KZ011-1562N081	8孔 8H	110.5	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws				

桡骨远端掌侧斜T型锁定板

Distal radius palmar oblique T-locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1572L030	3孔 3H	43	纯钛 Pure Titanium
KZ011-1572L040	4孔 4H	52	
KZ011-1572L050	5孔 5H	61	
KZ011-1572L060	6孔 6H	70	
KZ011-1572L070	7孔 7H	79	
KZ011-1572R030	3孔 3H	43	
KZ011-1572R040	4孔 4H	52	
KZ011-1572R050	5孔 5H	61	
KZ011-1572R060	6孔 6H	70	
KZ011-1572R070	7孔 7H	79	
备注：选用2.4/2.7锁定螺钉，2.4/2.7皮质骨螺钉 Remark:Select 2.4/2.7 locking screws, 2.4/2.7 cortical screws			

桡骨远端背侧径向柱锁定板

Distal radius dorsal radial column locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2202N050	5孔 5H	46	纯钛 Pure Titanium
KZ011-2202N060	6孔 6H	53	
KZ011-2202N070	7孔 7H	60	
KZ011-2202N080	8孔 8H	67	
备注：选用2.4/2.7万向锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 universal locking screws, 2.7 cortical screws.			

桡骨远端锁定板

Distal radius locking plate

桡骨远端关节外T型锁定板- I （头四孔）

Distal radius extra-articular T-locking plate-I (head four holes)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1542L030	3孔 3H	49	纯钛 Pure Titanium
KZ011-1542L040	4孔 4H	58	
KZ011-1542L050	5孔 5H	67	
KZ011-1542L060	6孔 6H	76	
KZ011-1542L070	7孔 7H	85	
KZ011-1542R030	3孔 3H	49	
KZ011-1542R040	4孔 4H	58	
KZ011-1542R050	5孔 5H	67	
KZ011-1542R060	6孔 6H	76	
KZ011-1542R070	7孔 7H	85	
备注：选用2.4/2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 locking screws, 2.7 cortical screws			

桡骨远端掌前侧锁定板- II

Distal radius palmar locking plate-II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1552L030	3孔 3H	49	纯钛 Pure Titanium
KZ011-1552L040	4孔 4H	58	
KZ011-1552L050	5孔 5H	67	
KZ011-1552L060	6孔 6H	76	
KZ011-1552L070	7孔 7H	85	
KZ011-1552R030	3孔 3H	49	
KZ011-1552R040	4孔 4H	58	
KZ011-1552R050	5孔 5H	67	
KZ011-1552R060	6孔 6H	76	
KZ011-1552R070	7孔 7H	85	
备注：选用2.4/2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 locking screws, 2.7 cortical screws			

桡骨远端掌侧双柱锁定板- I

Distal radius palmar double-column locking plate-I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1052L020	2孔 2H	42	纯钛 Pure Titanium
KZ011-1052L030	3孔 3H	51	
KZ011-1052L040	4孔 4H	63	
KZ011-1052L050	5孔 5H	72	
KZ011-1052R020	2孔 2H	42	
KZ011-1052R030	3孔 3H	51	
KZ011-1052R040	4孔 4H	63	
KZ011-1052R050	5孔 5H	72	
备注：选用2.4/2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 locking screws, 2.7 cortical screws			

桡骨远端锁定板

Distal radius locking plate

桡骨远端掌侧双柱锁定板- I

Distal radius palmar double-column locking plate-I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2022L020	2孔 2H	42	纯钛 Pure Titanium
KZ011-2022L030	3孔 3H	51	
KZ011-2022L040	4孔 4H	63	
KZ011-2022L050	5孔 5H	72	
KZ011-2022R020	2孔 2H	42	
KZ011-2022R030	3孔 3H	51	
KZ011-2022R040	4孔 4H	63	
KZ011-2022R050	5孔 5H	72	
备注：选用2.4/2.7万向锁定螺钉，2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 Universal Locking Screws, 2.7 Locking Screws, and 2.7 Cortical Screws.			

桡骨远端掌侧双柱锁定板- II

Distal radius palmar double-column locking plate-II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1042L020	2孔 2H	47	纯钛 Pure Titanium
KZ011-1042L030	3孔 3H	55	
KZ011-1042L040	4孔 4H	68	
KZ011-1042L050	5孔 5H	77	
KZ011-1042R020	2孔 2H	47	
KZ011-1042R030	3孔 3H	55	
KZ011-1042R040	4孔 4H	68	
KZ011-1042R050	5孔 5H	77	
备注：选用2.4/2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 locking screws, 2.7 cortical screws			

桡骨远端背侧中间柱锁定板

Dorsal intermediate column locking plate of distal radius

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2212L030	3孔 3H	41	纯钛 Pure Titanium
	KZ011-2212L040	4孔 4H	49	
	KZ011-2212R030	3孔 3H	41	
	KZ011-2212R040	4孔 4H	49	
	备注：选用2.4/2.7万向锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.4/2.7 universal locking screws, 2.7 cortical screws.			

桡骨远端锁定板

Distal radius locking plate

桡骨远端掌侧双柱锁定板- II

Distal radius palmar double-column locking plate-II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2012L020	2孔 2H	47	纯钛 Pure Titanium
KZ011-2012L030	3孔 3H	55	
KZ011-2012L040	4孔 4H	68	
KZ011-2012L050	5孔 5H	77	
KZ011-2012R020	2孔 2H	47	
KZ011-2012R030	3孔 3H	55	
KZ011-2012R040	4孔 4H	68	
KZ011-2012R050	5孔 5H	77	
备注：选用2.7锁定螺钉，2.4/2.7万向锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.4/2.7 polyaxial locking screws, and 2.7 cortical screws.			

桡骨远端内侧 I 型

Distal radial medial type I

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2162N030	3孔 3H	50.2	纯钛 Pure Titanium
KZ011-2162N040	4孔 4H	58.6	
KZ011-2162N050	5孔 5H	67	
KZ011-2162N070	7孔 7H	83.6	
KZ011-2162N090	9孔 9H	100.6	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

桡骨远端背侧IV型

Type IV distal radius dorsal

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2172N080	8孔 8H	132.5	纯钛 Pure Titanium
	KZ011-2172N100	10孔 10H	159.5	
	KZ011-2172N120	12孔 12H	186.5	
	备注：头部选用螺钉HC2.7，尾部选用螺钉HC3.5,HA3.5 Remark:Select HC2.7 for the head, HC3.5 and HA3.5 for the tail.			

桡骨远端内侧 II 型锁定板（左右式）

Type II locking plate on the medial aspect of the distal radius (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2182N040	4孔 4H	61.7	纯钛 Pure Titanium
	KZ011-2182N050	5孔 5H	71.7	
	KZ011-2182N060	6孔 6H	81.7	
	KZ011-2182N070	7孔 7H	91.7	
	KZ011-2182N080	8孔 8H	101.7	
	备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

桡骨远端锁定板

Distal radius locking plate

小T型斜锁定板（左右式）

Small T-shaped angled locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2192N030	3孔 3H	53	纯钛 Pure Titanium
KZ011-2192N040	4孔 4H	64	
KZ011-2192N050	5孔 5H	75	
KZ011-2192N060	6孔 6H	86	
KZ011-2192N070	7孔 7H	97	
备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark: Select 3.5 locking screws, 3.5 cortical screws.			

桡骨远端T型锁定板（头部4孔）

Distal radius T-locking plate (4 holes in the head)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2232N030	3孔 3H	50	纯钛 Pure Titanium
KZ011-2232N040	4孔 4H	62	
KZ011-2232N050	5孔 5H	74	
KZ011-2232N060	6孔 6H	86	
KZ011-2232N070	7孔 7H	98	
KZ011-2232N080	8孔 8H	110	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws, 3.5 cortical screws.			

桡骨远端钢板（左右式）

Distal radius plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2222L020	2孔 2H	56.2	纯钛 Pure Titanium
KZ011-2222L030	3孔 3H	64.7	
KZ011-2222L040	4孔 4H	73.2	
KZ011-2222L050	5孔 5H	81.7	
KZ011-2222L070	7孔 7H	96.7	
KZ011-2222L090	9孔 9H	115.7	
KZ011-2222R020	2孔 2H	56.2	
KZ011-2222R030	3孔 3H	64.7	
KZ011-2222R040	4孔 4H	73.2	
KZ011-2222R050	5孔 5H	81.7	
KZ011-2222R070	7孔 7H	96.7	
KZ011-2222R090	9孔 9H	115.7	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

桡骨远端锁定板

Distal radius locking plate

桡骨远端背侧锁定板Ⅶ（左右式）


Distal radius dorsal locking plate Ⅶ (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2242L030	3孔 3H	43	纯钛 Pure Titanium
KZ011-2242L040	4孔 4H	50	
KZ011-2242L050	5孔 5H	57	
KZ011-2242L060	6孔 6H	64	
KZ011-2242R030	3孔 3H	43	
KZ011-2242R040	4孔 4H	50	
KZ011-2242R050	5孔 5H	57	
KZ011-2242R060	6孔 6H	64	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉			
Remark:Select 2.7 locking screws, 2.7 cortical screws.			

小T型弯锁定板

Small T-shaped bend locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2252N030	3孔 3H	48	纯钛 Pure Titanium
KZ011-2252N040	4孔 4H	59	
KZ011-2252N050	5孔 5H	70	
KZ011-2252N060	6孔 6H	81	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉			
Remark:Select 3.5 locking screws, 3.5 cortical screws.			


2.4桡骨远端掌侧加长锁定板Ⅰ

2.4 Extended locking plate on the palmar side of the distal radiusⅠ

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.017.200	8孔 8H	100	纯钛 Pure Titanium
	101.017.204	10孔 10H	120	
	101.017.208	12孔 12H	140	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉			
	Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4桡骨远端掌侧关节外锁定板Ⅱ

2.4 Distal radius palmar external locking plateⅡ



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.025.200	3孔 3H	47	纯钛 Pure Titanium
101.025.201	3孔 3H	47	
101.025.202	4孔 4H	56	
101.025.203	4孔 4H	56	
101.025.204	5孔 5H	65	
101.025.205	5孔 5H	65	
备注：选用2.4锁定螺钉，2.4皮质骨螺钉 Remark:Select 2.4 locking screws,2.4 cortical screws			

桡骨远端锁定板

Distal radius locking plate

2.4桡骨远端掌侧加长锁定板Ⅱ

2.4 Extended locking plateⅡ on the palmar side of the distal radius



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.029.200	5孔 5H	95	纯钛 Pure Titanium
101.029.204	7孔 7H	125	
101.029.208	9孔 9H	154	
101.029.212	11孔 11H	184	
101.029.216	13孔 13H	212	
101.029.201	5孔 5H	95	
101.029.205	7孔 7H	125	
101.029.209	9孔 9H	154	
101.029.213	11孔 11H	184	
101.029.217	13孔 13H	212	
备注：选用2.4锁定螺钉，2.4皮质骨螺钉			
Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4桡骨头锁定板Ⅰ

2.4 Radial head locking plateⅠ

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	101.030.200	2孔 2H	32	纯钛 Pure Titanium
	101.030.201	3孔 3H	41	
	101.030.202	4孔 4H	50	
	备注：选用2.4锁定螺钉，2.4皮质骨螺钉			
	Remark:Select 2.4 locking screws,2.4 cortical screws			

2.4桡骨远端外侧锁定板

2.4 Lateral locking plate of distal radius



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.023.200	5孔 5H	48	纯钛 Pure Titanium
101.023.201	6孔 6H	57	
备注：选用2.7锁定螺钉，2.7皮质骨螺钉 Remark:Select 2.7 locking screws, 2.7 cortical screws.			

下肢锁定板Lower limb locking plate

股骨近端蛇型锁定板Proximal femoral snake locking plate



Product No.	Specification	Length	Material
产品编号	规格	长度(mm)	
KZ011-1102L040	4孔 4H	144	纯钛 Pure Titanium
KZ011-1102L050	5孔 5H	162	
KZ011-1102L060	6孔 6H	180	
KZ011-1102L070	7孔 7H	198	
KZ011-1102L080	8孔 8H	216	
KZ011-1102L090	9孔 9H	234	
KZ011-1102L100	10孔 10H	252	
KZ011-1102L120	12孔 12H	288	
KZ011-1102L140	14孔 14H	324	
KZ011-1102R040	4孔 4H	144	
KZ011-1102R050	5孔 5H	162	
KZ011-1102R060	6孔 6H	180	
KZ011-1102R070	7孔 7H	198	
KZ011-1102R080	8孔 8H	216	
KZ011-1102R090	9孔 9H	234	
KZ011-1102R100	10孔 10H	252	
KZ011-1102R120	12孔 12H	288	
KZ011-1102R140	14孔 14H	324	
备注：选用5.0锁定螺钉，6.5锁定接骨螺钉（空心），6.5锥形空心螺钉			
Remark:Select 5.0 locking screws, 6.5 locking cancellous screws (cannulated), 6.5 tapered cannulated screws.			

股骨近端蛇型锁定板Proximal femoral snake locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1962L040	4孔 4H	144	纯钛 Pure Titanium
KZ011-1962L050	5孔 5H	162	
KZ011-1962L060	6孔 6H	180	
KZ011-1962L070	7孔 7H	198	
KZ011-1962L080	8孔 8H	216	
KZ011-1962L90	9孔 9H	234	
KZ011-1962L100	10孔 10H	252	
KZ011-1962L120	12孔 12H	288	
KZ011-1962L140	14孔 14H	324	
KZ011-1962R040	4孔 4H	144	
KZ011-1962R050	5孔 5H	162	
KZ011-1962R060	6孔 6H	180	
KZ011-1962R070	7孔 7H	198	
KZ011-1962R080	8孔 8H	216	
KZ011-1962R090	9孔 9H	234	
KZ011-1962R100	10孔 10H	252	
KZ011-1962R120	12孔 12H	288	
KZ011-1962R140	14孔 14H	324	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉，6.5锁定接骨螺钉（空心），6.5锥形空心螺钉			
Remark:Select 5.0 locking screws,4.5 Cortical Screw,6.5 locking cancellous screws (cannulated), 6.5 tapered cannulated screws.			

下肢锁定板Lower limb locking plate

弯曲型锁定板Bending locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1222N080	8孔 8H	178.5	纯钛 Pure Titanium
KZ011-1222N100	10孔 10H	192	
KZ011-1222N110	11孔 11H	210	
KZ011-1222N120	12孔 12H	228	
KZ011-1222N130	13孔 13H	246	
KZ011-1222N140	14孔 14H	264	
KZ011-1222N150	15孔 15H	282	
KZ011-1222N160	16孔 16H	300	
KZ011-1222N180	18孔 18H	332	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

股骨有限接触锁定板Femoral Limited Contact Locking Plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1272N060	6孔 6H	116	纯钛 Pure Titanium
KZ011-1272N070	7孔 7H	134	
KZ011-1272N080	8孔 8H	152	
KZ011-1272N090	9孔 9H	170	
KZ011-1272N100	10孔 10H	188	
KZ011-1272N110	11孔 11H	206	
KZ011-1272N120	12孔 12H	224	
KZ011-1272N130	13孔 13H	242	
KZ011-1272N140	14孔 14H	260	
KZ011-1272N150	15孔 15H	278	
KZ011-1272N160	16孔 16H	296	
KZ011-1272N170	17孔 17H	314	
KZ011-1272N180	18孔 18H	332	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

股骨远端外侧锁定板Distal lateral femoral locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1192L050	5孔 5H	156	纯钛 Pure Titanium
KZ011-1192L070	7孔 7H	196	
KZ011-1192L090	9孔 9H	236	
KZ011-1192L110	11孔 11H	276	
KZ011-1192L130	13孔 13H	316	
KZ011-1192R050	5孔 5H	156	
KZ011-1192R070	7孔 7H	196	
KZ011-1192R090	9孔 9H	236	
KZ011-1192R110	11孔 11H	276	
KZ011-1192R130	13孔 13H	316	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉			
Remark:Select 5.0 locking screws,4.5 Cortical Screw			

下肢锁定板Lower limb locking plate

股骨远端解剖型锁定板Distal femoral anatomical locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1972L050	5孔 5H	142	纯钛 Pure Titanium
KZ011-1972L060	6孔 6H	160	
KZ011-1972L070	7孔 7H	178	
KZ011-1972L080	8孔 8H	196	
KZ011-1972L090	9孔 9H	214	
KZ011-1972L100	10孔 10H	232	
KZ011-1972L110	11孔 11H	250	
KZ011-1972L120	12孔 12H	268	
KZ011-1972L130	13孔 13H	286	
KZ011-1972R050	5孔 5H	142	
KZ011-1972R060	6孔 6H	160	
KZ011-1972R070	7孔 7H	178	
KZ011-1972R080	8孔 8H	196	
KZ011-1972R090	9孔 9H	214	
KZ011-1972R100	10孔 10H	232	
KZ011-1972R110	11孔 11H	250	
KZ011-1972R120	12孔 12H	268	
KZ011-1972R130	13孔 13H	286	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

股骨远端内侧截骨型锁定板Distal medial femoral osteotomy locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1632L040	4孔 4H	121	纯钛 Pure Titanium
	KZ011-1632R040	4孔 4H	121	
	备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

股骨远端外侧截骨型锁定板Distal lateral femoral osteotomy locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1642L050	5孔 5H	145	纯钛 Pure Titanium
KZ011-1642R050	5孔 5H	145	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

胫骨近端内侧截骨型锁定板Tibial medial osteotomy locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1612N041	4孔 4H	115	纯钛 Pure Titanium
	备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

下肢锁定板Lower limb locking plate

胫骨近端外侧截骨型锁定板Tibial lateral osteotomy locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1622L030	3孔 3H	102	纯钛 Pure Titanium
	KZ011-1622R030	3孔 3H	102	
	备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

胫骨近端外侧锁定板Lateral locking plate of proximal tibia



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1182L040	4孔 4H	123	纯钛 Pure Titanium
KZ011-1182L050	5孔 5H	140	
KZ011-1182L060	6孔 6H	160	
KZ011-1182L070	7孔 7H	180	
KZ011-1182L080	8孔 8H	200	
KZ011-1182L090	9孔 9H	220	
KZ011-1182L100	10孔 10H	240	
KZ011-1182L110	11孔 11H	260	
KZ011-1182L120	12孔 12H	280	
KZ011-1182L130	13孔 13H	300	
KZ011-1182R040	4孔 4H	123	
KZ011-1182R050	5孔 5H	140	
KZ011-1182R060	6孔 6H	160	
KZ011-1182R070	7孔 7H	180	
KZ011-1182R080	8孔 8H	200	
KZ011-1182R090	9孔 9H	220	
KZ011-1182R100	10孔 10H	240	
KZ011-1182R110	11孔 11H	260	
KZ011-1182R120	12孔 12H	280	
KZ011-1182R130	13孔 13H	300	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

胫骨近端L型锁定板- I L-shaped locking plate for proximal tibia - I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1122L040	4孔 4H	82	纯钛 Pure Titanium
KZ011-1122L060	6孔 6H	118	
KZ011-1122L080	8孔 8H	154	
KZ011-1122L100	10孔 10H	190	
KZ011-1122L120	12孔 12H	226	
KZ011-1122L140	14孔 14H	262	
KZ011-1122R040	4孔 4H	82	
KZ011-1122R060	6孔 6H	118	
KZ011-1122R080	8孔 8H	154	
KZ011-1122R100	10孔 10H	190	
KZ011-1122R120	12孔 12H	226	
KZ011-1122R140	14孔 14H	262	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

下肢锁定板Lower limb locking plate

胫骨近端L型锁定板（头七孔）L-shaped locking plate for proximal tibia (7-hole incision)

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1092L050	5孔 5H	92	纯钛 Pure Titanium
KZ011-1092L060	6孔 6H	104	
KZ011-1092L070	7孔 7H	116	
KZ011-1092L080	8孔 8H	128	
KZ011-1092L090	9孔 9H	140	
KZ011-1092L100	10孔 10H	152	
KZ011-1092L110	11孔 11H	164	
KZ011-1092L120	12孔 12H	176	
KZ011-1092L130	13孔 13H	188	
KZ011-1092L140	14孔 14H	200	
KZ011-1092L160	16孔 16H	224	
KZ011-1092R050	5孔 5H	92	
KZ011-1092R060	6孔 6H	104	
KZ011-1092R070	7孔 7H	116	
KZ011-1092R080	8孔 8H	128	
KZ011-1092R090	9孔 9H	140	
KZ011-1092R100	10孔 10H	152	
KZ011-1092R110	11孔 11H	164	
KZ011-1092R120	12孔 12H	176	
KZ011-1092R130	13孔 13H	188	
KZ011-1092R140	14孔 14H	200	
KZ011-1092R160	16孔 16H	224	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

胫骨近端T型锁定板T-locking plate for proximal tibia

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1372L040	4孔 4H	106	纯钛 Pure Titanium
KZ011-1372L060	6孔 6H	142	
KZ011-1372L080	8孔 8H	178	
KZ011-1372L100	10孔 10H	214	
KZ011-1372L120	12孔 12H	250	
KZ011-1372L140	14孔 14H	286	
KZ011-1372L160	16孔 16H	322	
KZ011-1372R040	4孔 4H	106	
KZ011-1372R060	6孔 6H	142	
KZ011-1372R080	8孔 8H	178	
KZ011-1372R100	10孔 10H	214	
KZ011-1372R120	12孔 12H	250	
KZ011-1372R140	14孔 14H	286	
KZ011-1372R160	16孔 16H	322	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

下肢锁定板Lower limb locking plate

胫骨近端匙内侧锁定板Tibial proximal key medial locking plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1142N020	2孔 2H	82	纯钛 Pure Titanium
KZ011-1142N030	3孔 3H	95	
KZ011-1142N040	4孔 4H	108	
KZ011-1142N050	5孔 5H	121	
KZ011-1142N060	6孔 6H	134	
KZ011-1142N070	7孔 7H	147	
KZ011-1142N080	8孔 8H	160	
KZ011-1142N090	9孔 9H	173	
KZ011-1142N100	10孔 10H	186	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

胫骨平台锁定板Tibial plateau locking plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1472L030	3孔 3H	103	纯钛 Pure Titanium
KZ011-1472L040	4孔 4H	121	
KZ011-1472L050	5孔 5H	139	
KZ011-1472L070	7孔 7H	175	
KZ011-1472L090	9孔 9H	211	
KZ011-1472L110	11孔 11H	247	
KZ011-1472L130	13孔 13H	283	
KZ011-1472R030	3孔 3H	103	
KZ011-1472R040	4孔 4H	121	
KZ011-1472R050	5孔 5H	139	
KZ011-1472R070	7孔 7H	175	
KZ011-1472R090	9孔 9H	211	
KZ011-1472R110	11孔 11H	247	
KZ011-1472R130	13孔 13H	283	
备注：选用5.0锁定螺钉，4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

不等厚干骺端锁定板Unequal thickness epiphyseal locking plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1232N030	3孔 3H	118	纯钛 Pure Titanium
KZ011-1232N040	4孔 4H	136	
KZ011-1232N050	5孔 5H	154	
KZ011-1232N060	6孔 6H	172	
KZ011-1232N070	7孔 7H	190	
KZ011-1232N080	8孔 8H	208	
KZ011-1232N090	9孔 9H	226	
KZ011-1232N110	11孔 11H	262	
KZ011-1232N130	13孔 13H	298	
KZ011-1232N150	15孔 15H	334	
备注：选用3.5/5.0锁定螺钉，3.5/4.5皮质骨螺钉 Remark:Select 3.5/5.0 locking screws,3.5/4.5 Cortical Screw			

下肢锁定板 Lower limb locking plate

胫骨直型锁定板-1 Tibial Straight Locking Plate-1


	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1252N060	6孔 6H	116	纯钛 Pure Titanium
	KZ011-1252N070	7孔 7H	134	
	KZ011-1252N080	8孔 8H	152	
	KZ011-1252N090	9孔 9H	170	
	KZ011-1252N100	10孔 10H	188	
	KZ011-1252N110	11孔 11H	206	
	KZ011-1252N120	12孔 12H	224	
	备注: 选用5.0锁定螺钉, 4.5皮质骨螺钉 Remark:Select 5.0 locking screws,4.5 Cortical Screw			

胫骨远端内侧锁定板-2 Distal Tibial Medial Locking Plate-2



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1132L040	4孔 4H	117	纯钛 Pure Titanium
KZ011-1132L050	5孔 5H	130.5	
KZ011-1132L060	6孔 6H	144	
KZ011-1132L070	7孔 7H	157.5	
KZ011-1132L080	8孔 8H	171	
KZ011-1132L090	9孔 9H	184.5	
KZ011-1132L100	10孔 10H	198	
KZ011-1132L110	11孔 11H	211.5	
KZ011-1132L120	12孔 12H	225	
KZ011-1132R040	4孔 4H	117	
KZ011-1132R050	5孔 5H	130.5	
KZ011-1132R060	6孔 6H	144	
KZ011-1132R070	7孔 7H	157.5	
KZ011-1132R080	8孔 8H	171	
KZ011-1132R090	9孔 9H	184.5	
KZ011-1132R100	10孔 10H	198	
KZ011-1132R110	11孔 11H	211.5	
KZ011-1132R120	12孔 12H	225	
备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

胫骨远端翼型锁定板 Distal tibial wing locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1012N050	5孔 5H	121	纯钛 Pure Titanium
	KZ011-1012N070	7孔 7H	147	
	KZ011-1012N090	9孔 9H	173	
	KZ011-1012N110	11孔 11H	199	
	备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

下肢锁定板 Lower limb locking plate

胫骨远端外侧L型锁定板 L-shaped locking plate on the distal lateral side of the tibia



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1162L050	5孔 5H	80	纯钛 Pure Titanium
KZ011-1162L070	7孔 7H	106	
KZ011-1162L090	9孔 9H	132	
KZ011-1162L110	11孔 11H	158	
KZ011-1162L130	13孔 13H	184	
KZ011-1162L150	15孔 15H	210	
KZ011-1162L170	17孔 17H	236	
KZ011-1162L190	19孔 19H	262	
KZ011-1162L210	21孔 21H	288	
KZ011-1162R050	5孔 5H	80	
KZ011-1162R070	7孔 7H	106	
KZ011-1162R090	9孔 9H	132	
KZ011-1162R110	11孔 11H	158	
KZ011-1162R130	13孔 13H	184	
KZ011-1162R150	15孔 15H	210	
KZ011-1162R170	17孔 17H	236	
KZ011-1162R190	19孔 19H	262	
KZ011-1162R210	21孔 21H	288	
备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

胫骨远端径内侧锁定板 I Medial locking plate of distal tibia I



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1152L040	4孔 4H	109	纯钛 Pure Titanium
KZ011-1152L050	5孔 5H	124	
KZ011-1152L060	6孔 6H	135	
KZ011-1152L070	7孔 7H	150	
KZ011-1152L080	8孔 8H	161	
KZ011-1152L090	9孔 9H	176	
KZ011-1152L100	10孔 10H	187	
KZ011-1152L110	11孔 11H	202	
KZ011-1152L120	12孔 12H	213	
KZ011-1152L130	13孔 13H	228	
KZ011-1152L140	14孔 14H	239	
KZ011-1152L150	15孔 15H	254	
KZ011-1152R040	4孔 4H	109	
KZ011-1152R050	5孔 5H	124	
KZ011-1152R060	6孔 6H	135	
KZ011-1152R070	7孔 7H	150	
KZ011-1152R080	8孔 8H	161	
KZ011-1152R090	9孔 9H	176	
KZ011-1152R100	10孔 10H	187	
KZ011-1152R110	11孔 11H	202	
KZ011-1152R120	12孔 12H	213	
KZ011-1152R130	13孔 13H	228	
KZ011-1152R140	14孔 14H	239	
KZ011-1152R150	15孔 15H	254	
备注: 选用3.5锁定螺钉, 3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

下肢锁定板 Lower limb locking plate

腓骨远端锁定板 Distal fibular locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1072L030	3孔 3H	73	纯钛 Pure Titanium
KZ011-1072L040	4孔 4H	86	
KZ011-1072L050	5孔 5H	99	
KZ011-1072L060	6孔 6H	112	
KZ011-1072L070	7孔 7H	125	
KZ011-1072L090	9孔 9H	151	
KZ011-1072L110	11孔 11H	177	
KZ011-1072R030	3孔 3H	73	
KZ011-1072R040	4孔 4H	86	
KZ011-1072R050	5孔 5H	99	
KZ011-1072R060	6孔 6H	112	
KZ011-1072R070	7孔 7H	125	
KZ011-1072R090	9孔 9H	151	
KZ011-1072R110	11孔 11H	177	
备注：选用2.7/3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 2.7/3.5 locking screws,3.5 Cortical Screw			

腓骨远端髌外侧锁定板 Lateral locking plate of distal fibular condyle



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-1702N030	3孔 3H	73	纯钛 Pure Titanium
KZ011-1702N040	4孔 4H	86	
KZ011-1702N050	5孔 5H	99	
KZ011-1702N060	6孔 6H	112	
KZ011-1702N070	7孔 7H	125	
KZ011-1702N080	8孔 8H	138	
KZ011-1702N090	9孔 9H	151	
KZ011-1702N100	10孔 10H	164	
KZ011-1702N110	11孔 11H	177	
备注：选用2.7/3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 2.7/3.5 locking screws,3.5 Cortical Screw			

胫骨平台外侧锁定板II型（左右式） Lateral locking plate of the tibial plateau type II (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2502L050	5孔 5H	87	纯钛 Pure Titanium
KZ011-2502L070	7孔 7H	116	
KZ011-2502L090	9孔 9H	144	
KZ011-2502L110	11孔 11H	172	
KZ011-2502L130	13孔 13H	200	
KZ011-2502L150	15孔 15H	228	
KZ011-2502R050	5孔 5H	87	
KZ011-2502R070	7孔 7H	116	
KZ011-2502R090	9孔 9H	144	
KZ011-2502R110	11孔 11H	172	
KZ011-2502R130	13孔 13H	200	
KZ011-2502R150	15孔 15H	228	
备注：选用螺钉HC3.5/HC4.0,HA3.5 Remark:Select screws HC3.5/4.0,HA3.5			

下肢锁定板 Lower limb locking plate

胫骨平台内侧L型锁定板（左右式） L-shaped Tibial plateau medial locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2512L040	4孔 4H	60	纯钛 Pure Titanium
KZ011-2512L060	6孔 6H	86	
KZ011-2512L080	8孔 8H	112	
KZ011-2512R040	4孔 4H	60	
KZ011-2512R060	6孔 6H	86	
KZ011-2512R080	8孔 8H	112	
备注：选用螺钉HC4.0 Remark:Select screws HC4.0			

胫骨平台内侧锁定板（左右式） Tibial plateau medial locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2522L030	3孔 3H	66.2	纯钛 Pure Titanium
KZ011-2522L040	4孔 4H	79.2	
KZ011-2522L050	5孔 5H	92.2	
KZ011-2522L060	6孔 6H	105.2	
KZ011-2522L070	7孔 7H	118.2	
KZ011-2522L080	8孔 8H	131.2	
KZ011-2522L090	9孔 9H	144.2	
KZ011-2522L100	10孔 10H	157.2	
KZ011-2522R030	3孔 3H	66.2	
KZ011-2522R040	4孔 4H	79.2	
KZ011-2522R050	5孔 5H	92.2	
KZ011-2522R060	6孔 6H	105.2	
KZ011-2522R070	7孔 7H	118.2	
KZ011-2522R080	8孔 8H	131.2	
KZ011-2522R090	9孔 9H	144.2	
KZ011-2522R100	10孔 10H	157.2	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			


胫骨平台内侧锁定板（左右式） Tibial plateau medial locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2532L050	5孔 5H	131	纯钛 Pure Titanium
KZ011-2532L070	7孔 7H	163	
KZ011-2532L090	9孔 9H	195	
KZ011-2532L110	11孔 11H	227	
KZ011-2532R050	5孔 5H	131	
KZ011-2532R070	7孔 7H	163	
KZ011-2532R090	9孔 9H	195	
KZ011-2532R110	11孔 11H	227	
备注：选用螺钉HC5.0,HA4.5 Remark:Select screws HC5.0,HA4.5			


下肢锁定板Lower limb locking plate

3.5/5.0 L型锁定板（左右式）3.5/5.0 L-shaped locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2542L050	5孔 5H	98.4	纯钛 Pure Titanium
KZ011-2542L060	6孔 6H	114.4	
KZ011-2542L070	7孔 7H	130.4	
KZ011-2542L080	8孔 8H	146.4	
KZ011-2542L090	9孔 9H	162.4	
KZ011-2542L100	10孔 10H	178.4	
KZ011-2542R050	5孔 5H	98.4	
KZ011-2542R060	6孔 6H	114.4	
KZ011-2542R070	7孔 7H	130.4	
KZ011-2542R080	8孔 8H	146.4	
KZ011-2542R090	9孔 9H	162.4	
KZ011-2542R100	10孔 10H	178.4	
备注：选用螺钉HC5.0,HA4.5或HC3.5.HA3.5			
Remark:Select screws HC5.0,HA4.5 or HC3.5.HA3.5			

重建型L型锁定板（左右式）Reconstruction type L-shaped locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2552L040	4孔 4H	88.8	纯钛 Pure Titanium
KZ011-2552L050	5孔 5H	104.8	
KZ011-2552L060	6孔 6H	120.8	
KZ011-2552L070	7孔 7H	136.8	
KZ011-2552L080	8孔 8H	157.8	
KZ011-2552L090	9孔 9H	168.8	
KZ011-2552L100	10孔 10H	184.8	
KZ011-2552R040	4孔 4H	88.8	
KZ011-2552R050	5孔 5H	104.8	
KZ011-2552R060	6孔 6H	120.8	
KZ011-2552R070	7孔 7H	136.8	
KZ011-2552R080	8孔 8H	157.8	
KZ011-2552R090	9孔 9H	168.8	
KZ011-2552R100	10孔 10H	184.8	
备注：选用螺钉HC5.0,HA4.5			
Remark:Select screws HC5.0,HA4.5			

T型锁定板T-shaped locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2562N040	4孔 4H	82	纯钛 Pure Titanium
	KZ011-2562N050	5孔 5H	98	
	KZ011-2562N060	6孔 6H	114	
	KZ011-2562N070	7孔 7H	130	
	KZ011-2562N080	8孔 8H	146	
	备注：选用螺钉HC5.0,HA4.5			
	Remark:Select screws HC5.0,HA4.5			

下肢锁定板Lower limb locking plate

胫骨平台后柱锁定板（左右式）Tibial plateau posterior column locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2572L030	3孔 3H	65	纯钛 Pure Titanium
KZ011-2572L050	5孔 5H	91	
KZ011-2572L070	7孔 7H	117	
KZ011-2572R030	3孔 3H	65	
KZ011-2572R050	5孔 5H	91	
KZ011-2572R070	7孔 7H	117	
备注：选用螺钉HC3.5,HC4.0 Remark:Select screws HC3.5,HC4.0			

重建型胫骨远端外侧锁定板（左右式）Reconstructed distal lateral tibial locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2582L050	5孔 5H	126.5	纯钛 Pure Titanium
KZ011-2582L070	7孔 7H	162.5	
KZ011-2582L090	9孔 9H	198.5	
KZ011-2582L110	11孔 11H	234.5	
KZ011-2582R050	5孔 5H	126.5	
KZ011-2582R070	7孔 7H	162.5	
KZ011-2582R090	9孔 9H	198.5	
KZ011-2582R110	11孔 11H	234.5	
备注：选用螺钉HC5.0			
Remark:Select screws HC5.0			

胫骨远端外侧锁定板（左右式）Lateral locking plate of distal tibia (left and right)

Product No. 产品编号	Specification 规格	Length 长度(mm)
KZ011-2592L050	5孔 5H	118
KZ011-2592L070	7孔 7H	154
KZ011-2592L090	9孔 9H	190
KZ011-2592L110	11孔 11H	226
KZ011-2592L130	13孔 13H	262
KZ011-2592R050	5孔 5H	118
KZ011-2592R070	7孔 7H	154
KZ011-2592R090	9孔 9H	190
KZ011-2592R110	11孔 11H	226
KZ011-2592R130	13孔 13H	262

备注：选用螺钉HC5.0.HA4.5
Remark:Select screws HC5.0.HA4.5

Material
材质

纯钛
Pure Titanium

斜L型锁定板（左右式）Slanted L-shaped locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2602L040	4孔 4H	70	纯钛 Pure Titanium
KZ011-2602R040	4孔 4H	70	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

下肢锁定板Lower limb locking plate

胫骨远端后髁锁定板（左右式）Distal tibial posterior condyle locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2612L030	3孔 3H	58.5	纯钛 Pure Titanium
KZ011-2612L040	4孔 4H	68.5	
KZ011-2612L050	5孔 5H	78.5	
KZ011-2612L060	6孔 6H	88.5	
KZ011-2612L070	7孔 7H	98.5	
KZ011-2612L080	8孔 8H	108.5	
KZ011-2612R030	3孔 3H	58.5	
KZ011-2612R040	4孔 4H	68.5	
KZ011-2612R050	5孔 5H	78.5	
KZ011-2612R060	6孔 6H	88.5	
KZ011-2612R070	7孔 7H	98.5	
KZ011-2612R080	8孔 8H	108.5	
备注：选用螺钉HC3.5或HC4.0,HA3.5或HB4.0 Remark:Select screws HC3.5 or HC4.0,HA3.5 or HB4.0			

微型T型锁定板Miniature T-type locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2622L040	4孔 4H	72.2	纯钛 Pure Titanium
	备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

胫骨远端内侧（重建型）锁定板Distal medial tibial locking plate (reconstruction type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2632L040	4孔 4H	89	纯钛 Pure Titanium
	KZ011-2632L060	6孔 6H	111	
	备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 Cortical Screw			

头部3.5/2.7 腓骨远端锁定板（左右式）Head 3.5/2.7 Distal fibular locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2642L040	4孔 4H	83.5	纯钛 Pure Titanium
KZ011-2642L050	5孔 5H	95.5	
KZ011-2642L060	6孔 6H	107.5	
KZ011-2642L070	7孔 7H	119.5	
KZ011-2642L080	8孔 8H	131.5	
KZ011-2642R040	4孔 4H	83.5	
KZ011-2642R050	5孔 5H	95.5	
KZ011-2642R060	6孔 6H	107.5	
KZ011-2642R070	7孔 7H	119.5	
KZ011-2642R080	8孔 8H	131.5	
备注：选用螺钉HC3.5.HA3.5.HC2.7			
Remark:Select screws HC3.5.HA3.5.HC2.7			

下肢锁定板Lower limb locking plate

腓骨远端I型锁定板（左右式）Distal fibular type I locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2652L030	3孔 3H	75.2	纯钛 Pure Titanium
KZ011-2652L040	4孔 4H	88.2	
KZ011-2652L050	5孔 5H	101.2	
KZ011-2652L060	6孔 6H	114.2	
KZ011-2652L070	7孔 7H	127.2	
KZ011-2652L080	8孔 8H	140.2	
KZ011-2652L090	9孔 9H	153.2	
KZ011-2652L100	10孔 10H	166.2	
KZ011-2652R030	3孔 3H	75.2	
KZ011-2652R040	4孔 4H	88.2	
KZ011-2652R050	5孔 5H	101.2	
KZ011-2652R060	6孔 6H	114.2	
KZ011-2652R070	7孔 7H	127.2	
KZ011-2652R080	8孔 8H	140.2	
KZ011-2652R090	9孔 9H	153.2	
KZ011-2652R110	10孔 10H	166.2	
备注：选用螺钉HC3.5.HA3.5.HC2.7			
Remark:Select screws HC3.5.HA3.5.HC2.7			

股骨近端髁锁定板Proximal femoral condyle locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2662L040	4孔 4H	145	纯钛 Pure Titanium
KZ011-2662L050	5孔 5H	164	
KZ011-2662L070	7孔 7H	202	
KZ011-2662L090	9孔 9H	240	
KZ011-2662L110	11孔 11H	278	
KZ011-2662L130	13孔 13H	316	
KZ011-2662R040	4孔 4H	145	
KZ011-2662R050	5孔 5H	164	
KZ011-2662R070	7孔 7H	202	
KZ011-2662R090	9孔 9H	240	
KZ011-2662R110	11孔 11H	278	
KZ011-2662R130	13孔 13H	316	
备注：选用螺钉HC5.0,HC6.5,HA4.5			
Remark:Select screws HC5.0,HC6.5,HA4.5			

130° 角度钢板130° angle steel plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2672N020	2孔 2H	/	纯钛
	KZ011-2672N030	3孔 3H	/	Pure Titanium

下肢锁定板

Lower limb locking plate

刀片钉

Blade nail



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ011-2683N075- KZ011-2683N130	10	75-130	钛合金 Titanium Alloy
备注：长度为5mm一档 Remark:Length increments of 5 mm			

拉力螺钉


Lag screw



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ011-2693N075- KZ011-2693N130	7.2	75-130	钛合金 Titanium Alloy
备注：长度为5mm一档 Remark:Length increments of 5 mm			

股骨近端锁定板（左右式）

Proximal femoral locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2702L050	5孔 5H	147	纯钛 Pure Titanium
KZ011-2702L070	7孔 7H	184	
KZ011-2702L090	9孔 9H	221	
KZ011-2702L110	11孔 11H	258	
KZ011-2702L130	13孔 13H	295	
KZ011-2702R050	5孔 5H	147	
KZ011-2702R070	7孔 7H	184	
KZ011-2702R090	9孔 9H	221	
KZ011-2702R110	11孔 11H	258	
KZ011-2702R130	13孔 13H	295	
备注：选用螺钉HC5.0,HC6.5(实心或空心),HA4.5 Remark:Select screws HC5.0,HC6.5(Solid or Hollow),HA4.5			

股骨近端髁锁定板（左右式）

Proximal femoral condyle locking plate (left and right type)




Product No. 产品编号	Specification 规格	Thickness 厚度(mm)	Material 材质
KZ011-2722L040	/	5.6	纯钛 Pure Titanium
KZ011-2722R040	/	5.6	
备注：选用7.3空心锁定螺钉 Remark:Select 7.3 mm cannulated locking screws			

下肢锁定板

Lower limb locking plate

135° 角度板（左右式）

135° Angle Plate (Left and Right Type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2712L030	3孔 3H	158.2	纯钛 Pure Titanium
KZ011-2712L040	4孔 4H	183.2	
KZ011-2712L050	5孔 5H	208.2	
KZ011-2712L060	6孔 6H	233.2	
KZ011-2712L070	7孔 7H	258.2	
KZ011-2712L080	8孔 8H	283.2	
KZ011-2712L090	9孔 9H	308.2	
KZ011-2712L100	10孔 10H	333.2	
KZ011-2712L110	11孔 11H	358.2	
KZ011-2712L120	12孔 12H	383.2	
KZ011-2712L130	13孔 13H	408.2	
KZ011-2712L140	14孔 14H	433.2	
KZ011-2712L150	15孔 15H	456.2	
KZ011-2712L160	16孔 16H	483.2	
KZ011-2712R030	3孔 3H	158.2	
KZ011-2712R040	4孔 4H	183.2	
KZ011-2712R050	5孔 5H	208.2	
KZ011-2712R060	6孔 6H	233.2	
KZ011-2712R070	7孔 7H	258.2	
KZ011-2712R080	8孔 8H	283.2	
KZ011-2712R090	9孔 9H	308.2	
KZ011-2712R100	10孔 10H	333.2	
KZ011-2712R110	11孔 11H	358.2	
KZ011-2712R120	12孔 12H	383.2	
KZ011-2712R130	13孔 13H	408.2	
KZ011-2712R140	14孔 14H	433.2	
KZ011-2712R150	15孔 15H	456.2	
KZ011-2712R160	16孔 16H	483.2	
备注：选用螺钉HC5.0,HA4.5 Remark:Select screws HC5.0,HA4.5			

股骨近端髁锁定板 II 型（左右式）

Type II proximal femoral condyle locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2732L040	4孔 4H	137.4	纯钛 Pure Titanium
KZ011-2732L050	5孔 5H	156.4	
KZ011-2732L070	7孔 7H	194.4	
KZ011-2732L090	9孔 9H	232.4	
KZ011-2732L110	11孔 11H	270.4	
KZ011-2732L130	13孔 13H	308.4	
KZ011-2732R040	4孔 4H	137.4	
KZ011-2732R050	5孔 5H	156.4	
KZ011-2732R070	7孔 7H	194.4	
KZ011-2732R090	9孔 9H	232.4	
KZ011-2732R110	11孔 11H	270.4	
KZ011-2732R130	13孔 13H	308.4	
备注：选用螺钉HC5.0,HC6.5,HA4.5,HB6.5 Remark:Select screws HC5.0,HC6.5,HA4.5,HB6.5			

下肢锁定板

Lower limb locking plate

股骨近端髁锁定板III型（左右式）Type III proximal femoral condyle locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2742L020	2孔 2H	123	纯钛 Pure Titanium
KZ011-2742L050	5孔 5H	159	
KZ011-2742L060	6孔 6H	195	
KZ011-2742L080	8孔 8H	231	
KZ011-2742L100	10孔 10H	267	
KZ011-2742L120	12孔 12H	303	
KZ011-2742R020	2孔 2H	123	
KZ011-2742R050	5孔 5H	159	
KZ011-2742R060	6孔 6H	195	
KZ011-2742R080	8孔 8H	231	
KZ011-2742R100	10孔 10H	267	
KZ011-2742R120	12孔 12H	303	
备注：选用螺钉HC5.0,HC6.5,HA4.5			
Remark:Select screws HC5.0,HC6.5,HA4.5			

股骨近端髁锁定板 I 型（左右式）Type I proximal femoral condyle locking plate (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2752N050	5孔 5H	156.4	纯钛 Pure Titanium
	KZ011-2752N070	7孔 7H	194.4	
	KZ011-2752N090	9孔 9H	232.4	
	KZ011-2752N110	11孔 11H	270.4	
备注：选用螺钉HC5.0,HA4.5,HB6.5				
Remark:Select screws HC5.0,HA4.5,HB6.5				

股骨远端外侧髁锁定板（左右式）Distal femoral lateral condyle locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2762L030	3孔 3H	93	纯钛 Pure Titanium
KZ011-2762L050	5孔 5H	129.5	
KZ011-2762L070	7孔 7H	161.5	
KZ011-2762L090	9孔 9H	193.5	
KZ011-2762R030	3孔 3H	93	
KZ011-2762R050	5孔 5H	129.5	
KZ011-2762R070	7孔 7H	161.5	
KZ011-2762R090	9孔 9H	193.5	
备注：选用螺钉HC5.0,HA4.5 Remark:Select screws HC5.0,HA4.5			

下肢锁定板

Lower limb locking plate

股骨远端内侧髁锁定板（左右式）Distal femoral medial condyle locking plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2772L030	3孔 3H	97	纯钛 Pure Titanium
KZ011-2772L050	5孔 5H	125	
KZ011-2772L070	7孔 7H	157	
KZ011-2772L090	9孔 9H	188.5	
KZ011-2772R030	3孔 3H	97	
KZ011-2772R050	5孔 5H	125	
KZ011-2772R070	7孔 7H	157	
KZ011-2772R090	9孔 9H	188.5	
备注：选用螺钉HC5.0,HA4.5 Remark:Select screws HC5.0,HA4.5			

角度锁定板Angle locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2782L030	3孔 3H	113	纯钛 Pure Titanium
	KZ011-2782L040	4孔 4H	131	
	KZ011-2782L050	5孔 5H	149	
	KZ011-2782L060	6孔 6H	167	
	KZ011-2782L070	7孔 7H	185	
	KZ011-2782L080	8孔 8H	203	
	KZ011-2782L090	9孔 9H	221	
	KZ011-2782L100	10孔 10H	239	
	KZ011-2782L110	11孔 11H	257	
	KZ011-2782L120	12孔 12H	275	
	KZ011-2782R030	3孔 3H	113	
	KZ011-2782R040	4孔 4H	131	
	KZ011-2782R050	5孔 5H	149	
	KZ011-2782R060	6孔 6H	167	
	KZ011-2782R070	7孔 7H	185	
	KZ011-2782R080	8孔 8H	203	
	KZ011-2782R090	9孔 9H	221	
	KZ011-2782R100	10孔 10H	239	
	KZ011-2782R110	11孔 11H	257	
	KZ011-2782R120	12孔 12H	275	
备注：选用螺钉HC5.0,HA4.5,HB6.5				
Remark:Select screws HC5.0,HA4.5,HB6.5				

3.5腓骨远端锁定接骨板 II3.5 Distal fibular locking plate II



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.051.204	5孔 5H	99	纯钛 Pure Titanium
101.051.208	7孔 7H	125	
101.051.212	9孔 9H	151	
101.051.216	11孔 11H	177	
101.051.205	5孔 5H	99	
101.051.209	7孔 7H	125	
101.051.213	9孔 9H	151	
101.051.217	11孔 11H	177	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉 Remark:Select 3.5 locking screws,3.5 cortical screws			

下肢锁定板Lower limb locking plate

3.5胫骨近端外侧锁定接骨板Ⅰ3.5 Lateral locking plate of the proximal tibiaⅠ



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.053.200	4孔 4H	81	纯钛 Pure Titanium
101.053.204	6孔 6H	107	
101.053.208	8孔 8H	133	
101.053.212	10孔 10H	159	
101.053.216	12孔 12H	185	
101.053.220	14孔 14H	211	
101.053.201	4孔 4H	81	
101.053.205	6孔 6H	107	
101.053.209	8孔 8H	133	
101.053.213	10孔 10H	159	
101.053.217	12孔 12H	185	
101.053.221	14孔 14H	211	
备注：选用3.5锁定螺钉，3.5皮质骨螺钉			
Remark:Select 3.5 locking screws,3.5 cortical screws			


5.0胫骨近端外侧锁定接骨板Ⅲ5.0 Lateral locking plate of the proximal tibiaⅢ



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
101.062.202	5孔 5H	140	纯钛 Pure Titanium
101.062.206	7孔 7H	180	
101.062.210	9孔 9H	220	
101.062.214	11孔 11H	260	
101.062.203	5孔 5H	140	
101.062.207	7孔 7H	180	
101.062.211	9孔 9H	220	
101.062.215	11孔 11H	260	
备注：选用5.0锁定螺钉，5.0皮质骨螺钉 Remark:Select 5.0 locking screws,5.0 cortical screws			

微型锁定板Minisize locking plate

钛网Titanium mesh

	Product No. 产品编号	Length 长度(mm)	Material 材质
	KZ011-2792N080	80	纯钛 Pure Titanium
	KZ011-2792N090	90	
	KZ011-2792N100	100	
	KZ011-2792N120	120	
	KZ011-2792N140	140	
	KZ011-2792N150	150	
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw		

Z型钢板（左右式）Z-shaped plate (left and right type)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2802L040	4孔 4H	4.5	纯钛 Pure Titanium
	KZ011-2802R040	4孔 4H	4.5	
	备注：选用螺钉HA2.0（十字）			
	Remark:Select HA2.0 (Phillips) screw			

弓Ⅰ型钢板Bow-shaped plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2812N040	4孔 4H	23.7	纯 钛 Pure Titanium
KZ011-2812N050	5孔 5H	28.2	
KZ011-2812N060	6孔 6H	32.7	
KZ011-2812N070	7孔 7H	37.2	
KZ011-2812N080	8孔 8H	41.2	
KZ011-2812N090	9孔 9H	46.2	
KZ011-2812N100	10孔 10H	50.7	
KZ011-2812N110	11孔 11H	55.2	
KZ011-2812N120	12孔 12H	59.7	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

弓Ⅱ型钢板BowⅡtype plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2822N040	4孔 4H	22.2	纯 钛 Pure Titanium
	KZ011-2822N050	5孔 5H	27.9	
	KZ011-2822N060	6孔 6H	33.6	
	KZ011-2822N070	7孔 7H	39	
	KZ011-2822N080	8孔 8H	44.2	
	KZ011-2822N090	9孔 9H	49.2	
	KZ011-2822N100	10孔 10H	53.9	
	KZ011-2822N110	11孔 11H	58.3	
	KZ011-2822N120	12孔 12H	62.4	
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

微型锁定板Minisize locking plate

直型钢板（Ⅰ型）Straight plate (Type I)




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2832N040	4孔 4H	45.5	纯 钛 Pure Titanium
KZ011-2832N041	4孔 4H	45.5	
KZ011-2832N042	4孔 4H	45.5	
KZ011-2832N043	4孔 4H	45.5	
KZ011-2832N050	5孔 5H	45.5	
KZ011-2832N060	6孔 6H	45.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

直型钢板（Ⅱ型）Straight plate (Type II)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2842N040	4孔 4H	22.5	纯 钛 Pure Titanium
KZ011-2842N060	6孔 6H	34.5	
KZ011-2842N080	8孔 8H	46.5	
KZ011-2842N100	10孔 10H	58.5	
KZ011-2842N120	12孔 12H	70.5	
KZ011-2842N140	14孔 14H	82.5	
KZ011-2842N160	16孔 16H	94.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

斜L型钢板（左右式）Slanted L-shaped plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2852L050	5孔 5H	4.5	纯 钛 Pure Titanium
KZ011-2852R050	5孔 5H	4.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			


弯目型钢板CM Type plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2862N080	8孔 8H	27.1	纯 钛 Pure Titanium
KZ011-2862N120	12孔 12H	39.6	
KZ011-2862N160	16孔 16H	52.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			


微型锁定板Minisize locking plate

矩形圆孔型钢板ZM Type plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2872N040	4孔 4H	10.5	纯 钛 Pure Titanium
KZ011-2872N080	8孔 8H	22.5	
KZ011-2872N120	12孔 12H	34.5	
KZ011-2872N160	16孔 16H	46.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

正目型钢板YM Type plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2882N040	4孔 4H	10.5	纯 钛 Pure Titanium
KZ011-2882N080	8孔 8H	22.5	
KZ011-2882N120	12孔 12H	34.5	
KZ011-2882N160	16孔 16H	46.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

L型钢板Ⅰ型（左右式）L-shaped Plates Type I (Left and Right Types)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2892L040	4孔 4H	4.5	纯 钛 Pure Titanium
KZ011-2892L050	5孔 5H	4.5	
KZ011-2892L060	6孔 6H	4.5	
KZ011-2892R040	4孔 4H	4.5	
KZ011-2892R050	5孔 5H	4.5	
KZ011-2892R060	6孔 6H	4.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

L型钢板Ⅱ型（左右式）L-shaped Plates Type II (Left and Right Types)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-2902L040	4孔 4H	4.5	纯 钛 Pure Titanium
KZ011-2902L050	5孔 5H	4.5	
KZ011-2902L060	6孔 6H	4.5	
KZ011-2902R040	4孔 4H	4.5	
KZ011-2902R050	5孔 5H	4.5	
KZ011-2902R060	6孔 6H	4.5	
备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

微型锁定板Minisize locking plate

L型110° 钢板（左右式）L-shaped 110° Plates (Left and Right Types)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2912L040	4孔 4H	4.5	纯钛 Pure Titanium
	KZ011-2912L050	5孔 5H	4.5	
	KZ011-2912L060	6孔 6H	4.5	
	KZ011-2912R040	4孔 4H	4.5	
	KZ011-2912R050	5孔 5H	4.5	
	KZ011-2912R060	6孔 6H	4.5	
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

X型钢板X-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2922N060	6孔 6H	21.9	纯钛 Pure Titanium
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

叉型钢板Y-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2932N040	4孔 4H	14.7	纯钛 Pure Titanium
	KZ011-2932N050	5孔 5H	20.7	
	KZ011-2932N060	6孔 6H	26.7	
	KZ011-2932N070	7孔 7H	32.7	
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

T型钢板T-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2942N050	5孔 5H	23.5	纯钛 Pure Titanium
	KZ011-2942N060	6孔 6H	29.5	
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

工型钢板“工”-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2952N050	5孔 5H	18.3	纯钛 Pure Titanium
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

微型锁定板Minisize locking plate

H型钢板H-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2962N110	11孔 11H	28.5	纯钛 Pure Titanium
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

DY型钢板DY-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2972N120	12孔 12H	25	纯钛 Pure Titanium
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

DF型钢板DF-Shaped Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-2982N120	12孔 12H	48.7	纯 钛 Pure Titanium
	备注：选用螺钉HA2.0（十字） Remark:Select HA2.0 (Phillips) screw			

足踝锁定板Ankle locking plate

趾骨X型锁定板Toe X-shaped locking plate


	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1432N041	4孔 4H	18	纯钛 Pure Titanium
	KZ011-1432N042	4孔 4H	24	
	KZ011-1432N043	4孔 4H	30	
	备注：选用2.7锁定螺钉 Remark:Select 2.7 locking screws			

趾骨I型锁定板Toe I-shaped locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-1862N020	2孔 2H	17	纯钛 Pure Titanium
	备注：选用2.7锁定螺钉 Remark:Select 2.7 locking screws			


足踝锁定板 Ankle locking plate

跟骨锁定板 Calcaneal Locking Plate




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3202L080	8孔 8H	60	纯钛 Pure Titanium
KZ011-3202R080	8孔 8H	60	
备注: 选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

跟骨平台锁定板 Calcaneal plate locking plate




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3242L151	15孔 15H	64	纯钛 Pure Titanium
KZ011-3242L152	15孔 15H	69	
KZ011-3242L153	15孔 15H	76	
KZ011-3242R151	15孔 15H	64	
KZ011-3242R152	15孔 15H	69	
KZ011-3242R153	15孔 15H	76	
备注: 选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

跟骨长叶型锁定板 Calcaneal long-leaf locking plate (left and right)




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3212L150	15孔 15H	89	纯钛 Pure Titanium
KZ011-3212R150	15孔 15H	89	
备注: 选用3.5锁定螺钉 Remark:Select 3.5 locking screws			

跟骨环型锁定板 Calcaneal ring locking plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3262L120	12孔 12H	60	纯钛 Pure Titanium
KZ011-3262L121	12孔 12H	68	
KZ011-3262R120	12孔 12H	60	
KZ011-3262R121	12孔 12H	68	
备注: 选用3.5锁定螺钉 Remark:Select 3.5 locking screws			


跟骨短叶型锁定板 Calcaneal short-leaf locking plate



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3222L110	11孔	60	纯钛 Pure Titanium
KZ011-3222R110	11孔	60	
备注: 选用3.5锁定螺钉 Remark:Select 3.5 locking screws			


足踝锁定板 Ankle locking plate

舟骨锁定板 Navicular locking plate




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3232N110	11孔 11H	57	纯钛 Pure Titanium
备注: 选用2.7锁定螺钉 Remark:Select 2.7 locking screws			

跟骨V型锁定板（左右式） Calcaneal Type V Locking Plate (Left and Right)




Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3272L120	大Large	67.2	纯钛 Pure Titanium
KZ011-3272R120	大Large	67.2	
KZ011-3272L121	中Medium	58	
KZ011-3272R121	中Medium	58	
KZ011-3272L122	小Small	51	
KZ011-3272R122	小Small	51	
备注: 选用螺钉HC4.0 Remark:Select HC4.0 screws			

跟骨钢板V型解剖板（左右式） Calcaneal plate V-shaped anatomical plate (left and right type)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3282L110	大Large	72.6	纯钛 Pure Titanium
KZ011-3282R110	大Large	72.6	
KZ011-3282L111	中Medium	67.4	
KZ011-3282R111	中Medium	67.4	
KZ011-3282L112	小Small	60.3	
KZ011-3282R112	小Small	60.3	
备注: 选用螺钉HC3.5 Remark:Select HC3.5 screws			


跟骨IV型锁定板（左右式） Calcaneal Type IV Locking Plate (Left and Right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3292L120	大Large	58.8	纯钛 Pure Titanium
KZ011-3292R120	大Large	58.8	
KZ011-3292L121	小Small	48.8	
KZ011-32R121	小Small	48.8	
备注: 选用螺钉HC4.0 Remark:Select HC4.0 screws			

足踝锁定板Ankle locking plate

跟骨II型锁定板（左右式）Calcaneal Type II Locking Plate (Left and Right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3302L130	长3孔 Long 3H	79.5	纯钛 Pure Titanium
	KZ011-3302R130	长3孔 Long 3H	79.5	
	KZ011-3302L131	长4孔 Long 4H	89	
	KZ011-3302R131	长4孔 Long 4H	89	
	KZ011-3302L132	短3孔 Short 3H	59.3	
	KZ011-3302R132	短3孔 Short 3H	59.3	
	KZ011-3302L133	短4孔 Short 4H	69	
	KZ011-3302R133	短4孔 Short 4H	69	
备注：选用螺钉HC4.0 Remark:Select HC4.0 screws				

跟骨 I 型锁定板Calcaneal Type I Locking Plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3312L110	大 Large	63.3	纯钛 Pure Titanium
	KZ011-3312R110	大 Large	63.3	
	KZ011-3312L111	小 Small	49.5	
	KZ011-3312R111	小 Small	49.5	
备注：选用螺钉HC4.0 Remark:Select HC4.0 screws				

骰骨锁定板（左右式）Cuboid Locking Plate (Left and Right Types)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3322L090	9孔 9H	27.5	纯钛 Pure Titanium
	KZ011-3322R090	9孔 9H	27.5	
备注：选用螺钉 HC2.9 Remark:Select HC2.9screws				

舟骨锁定板（左右式）Navicular Bone locking plate (Left and Right Types)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3332L090	9孔 9H	50	纯钛 Pure Titanium
	KZ011-3332R090	9孔 9H	50	
备注：选用螺钉 HC2.9 Remark:Select HC2.9screws				

足踝锁定板Ankle locking plate

跟骨前突板（左右式）Anterior calcaneal plate (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3342L050	5孔 5H	50.7	纯钛 Pure Titanium
	KZ011-3342R050	5孔 5H	50.7	
备注：选用螺钉 HC3.5 Remark:Select HC3.5screws				

跖骨T型板Metatarsal bone T-Type locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3352N020	2孔 2H	27.4	纯钛 Pure Titanium
	KZ011-3352N030	3孔 3H	33.9	
	KZ011-3352N040	4孔 4H	40.4	
	KZ011-3352N050	5孔 5H	47	
	KZ011-3352N060	6孔 6H	53.7	
备注：选用3.5锁定螺钉 Remark:Select 3.5 locking screws				

跖趾融合锁定板（左右式）Metatarsophalangeal fusion locking plate (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3362L070	7孔 7H	56	纯钛 Pure Titanium
	KZ011-3362R070	7孔 7H	56	
备注：选用螺钉 HC2.9、HA2.7(小帽) Remark:Select HC2.9 screws、HA2.7(small nut) Screws				

距骨锁定板Talus locking plate

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3372N060	6孔 6H	38.5	纯钛 Pure Titanium
	备注：选用螺钉 HC2.7 Remark:Select HC2.7 screws			

跟骨结合型锁定板（左右式）Calcaneal-integrated locking plate (left and right)

	Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
	KZ011-3382L080	8孔 8H	56.9	纯钛 Pure Titanium
	KZ011-3382R080	8孔 8H	56.9	
	KZ011-3382L090	9孔 9H	68.7	
	KZ011-3382R090	9孔 9H	68.7	
备注：选用螺钉 HC3.5 Remark:Select HC3.5screws				

足踝锁定板 Ankle locking plate

跟骨后粗隆板（左右式） Posterior calcaneal plate (left and right)

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3392L060	6孔 6H	49.2	纯钛 Pure Titanium
KZ011-3392R060	6孔 6H	49.2	
KZ011-3392L070	7孔 7H	58.7	
KZ011-3392R070	7孔 7H	58.7	
备注：选用螺钉 HC3.5 Remark:Select HC3.5 screws			

距骨板 Talus plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3402N050	5孔 5H	43	纯钛 Pure Titanium
备注：选用螺钉 HC2.9 Remark:Select HC2.9 screws			

距骨近端解剖板（左右式） Proximal talus anatomical plate (left and right)



Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3412L030	3孔 3H	38	纯钛 Pure Titanium
KZ011-3412R030	3孔 3H	38	
KZ011-3412L040	4孔 4H	45	
KZ011-3412R040	4孔 4H	45	
KZ011-3412L050	5孔 5H	51.5	
KZ011-3412R050	5孔 5H	51.5	
KZ011-3412L060	6孔 6H	58	
KZ011-3412R060	6孔 6H	58	
KZ011-3412L070	7孔 7H	64.5	
KZ011-3412R070	7孔 7H	64.5	
KZ011-3412L080	8孔 8H	71	
KZ011-3412R080	8孔 8H	71	
备注：选用螺钉 HC3.5 Remark:Select HC3.5 screws			

跖楔融合板 Plantar wedge fusion plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3422N080	8孔 8H	60	纯钛 Pure Titanium
备注：选用螺钉 HC2.9、HA2.7(小帽) Remark:Select HC 2.9screws、HA2.7(small nut) Screws			

X形板 X-shaped plate

Product No. 产品编号	Specification 规格	Length 长度(mm)	Material 材质
KZ011-3432N080	4孔 4H	60	纯钛 Pure Titanium
备注：选用螺钉 HC3.5 Remark:Select HC3.5 screws			

螺钉系列 Screw series

2.4锁定螺钉（万向） 2.4 Locking screws (universal)


Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1023N060	2.4	6	钛合金 Titanium Alloy
KZ012-1023N080	2.4	8	
KZ012-1023N100	2.4	10	
KZ012-1023N120	2.4	12	
KZ012-1023N140	2.4	14	
KZ012-1023N160	2.4	16	
KZ012-1023N180	2.4	18	
KZ012-1023N200	2.4	20	
KZ012-1023N220	2.4	22	
KZ012-1023N240	2.4	24	
KZ012-1023N260	2.4	26	
KZ012-1023N280	2.4	28	
KZ012-1023N300	2.4	30	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

2.7锁定螺钉（万向） 2.7 Locking screws (universal)

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1013N060	2.7	6	钛合金 Titanium Alloy
KZ012-1013N080	2.7	8	
KZ012-1013N100	2.7	10	
KZ012-1013N120	2.7	12	
KZ012-1013N140	2.7	14	
KZ012-1013N160	2.7	16	
KZ012-1013N180	2.7	18	
KZ012-1013N200	2.7	20	
KZ012-1013N220	2.7	22	
KZ012-1013N240	2.7	24	
KZ012-1013N260	2.7	26	
KZ012-1013N280	2.7	28	
KZ012-1013N300	2.7	30	
KZ012-1013N320	2.7	32	
KZ012-1013N340	2.7	34	
KZ012-1013N360	2.7	36	
KZ012-1013N380	2.7	38	
KZ012-1013N400	2.7	40	
KZ012-1013N460	2.7	46	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列 Screw series

1.5锁定螺钉 1.5 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1033N060	1.5	6	钛合金 Titanium Alloy
KZ012-1033N070	1.5	7	
KZ012-1033N080	1.5	8	
KZ012-1033N090	1.5	9	
KZ012-1033N100	1.5	10	
KZ012-1033N110	1.5	11	
KZ012-1033N120	1.5	12	
KZ012-1033N130	1.5	13	
KZ012-1033N140	1.5	14	
KZ012-1033N150	1.5	15	
KZ012-1033N160	1.5	16	
KZ012-1033N180	1.5	18	
KZ012-1033N200	1.5	20	
KZ012-1033N220	1.5	22	
KZ012-1033N240	1.5	24	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

2.0锁定螺钉 2.0 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1043N060	2.0	6	钛合金 Titanium Alloy
KZ012-1043N080	2.0	8	
KZ012-1043N100	2.0	10	
KZ012-1043N120	2.0	12	
KZ012-1043N140	2.0	14	
KZ012-1043N160	2.0	16	
KZ012-1043N180	2.0	18	
KZ012-1043N200	2.0	20	
KZ012-1043N220	2.0	22	
KZ012-1043N240	2.0	24	
KZ012-1043N260	2.0	26	
KZ012-1043N280	2.0	28	
KZ012-1043N300	2.0	30	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列 Screw series

2.4锁定螺钉 2.4 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1093N060	2.4	6	钛合金 Titanium Alloy
KZ012-1093N080	2.4	8	
KZ012-1093N100	2.4	10	
KZ012-1093N120	2.4	12	
KZ012-1093N140	2.4	14	
KZ012-1093N160	2.4	16	
KZ012-1093N180	2.4	18	
KZ012-1093N200	2.4	20	
KZ012-1093N220	2.4	22	
KZ012-1093N240	2.4	24	
KZ012-1093N260	2.4	26	
KZ012-1093N280	2.4	28	
KZ012-1093N300	2.4	30	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

2.7锁定螺钉 2.7 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1103N060	2.7	6	钛合金 Titanium Alloy
KZ012-1103N080	2.7	8	
KZ012-1103N100	2.7	10	
KZ012-1103N120	2.7	12	
KZ012-1103N140	2.7	14	
KZ012-1103N160	2.7	16	
KZ012-1103N180	2.7	18	
KZ012-1103N200	2.7	20	
KZ012-1103N220	2.7	22	
KZ012-1103N240	2.7	24	
KZ012-1103N260	2.7	26	
KZ012-1103N280	2.7	28	
KZ012-1103N300	2.7	30	
KZ012-1103N320	2.7	32	
KZ012-1103N340	2.7	34	
KZ012-1103N360	2.7	36	
KZ012-1103N380	2.7	38	
KZ012-1103N400	2.7	40	
KZ012-1103N420	2.7	42	
KZ012-1103N440	2.7	44	
KZ012-1103N460	2.7	46	
KZ012-1103N480	2.7	48	
KZ012-1103N500	2.7	50	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列Screw series

2.9锁定螺钉2.9 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1193N060	2.9	6	钛合金 Titanium Alloy
KZ012-1193N080	2.9	8	
KZ012-1193N100	2.9	10	
KZ012-1193N120	2.9	12	
KZ012-1193N140	2.9	14	
KZ012-1193N160	2.9	16	
KZ012-1193N180	2.9	18	
KZ012-1193N200	2.9	20	
KZ012-1193N220	2.9	22	
KZ012-1193N240	2.9	24	
KZ012-1193N260	2.9	26	
KZ012-1193N280	2.9	28	
KZ012-1193N300	2.9	30	
KZ012-1193N320	2.9	32	
KZ012-1193N340	2.9	34	
KZ012-1193N360	2.9	36	
KZ012-1193N380	2.9	38	
KZ012-1193N400	2.9	40	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

3.5锁定螺钉3.5 Locking screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1113N100	3.5	10	钛合金 Titanium Alloy
KZ012-1113N120	3.5	12	
KZ012-1113N140	3.5	14	
KZ012-1113N160	3.5	16	
KZ012-1113N180	3.5	18	
KZ012-1113N200	3.5	20	
KZ012-1113N220	3.5	22	
KZ012-1113N240	3.5	24	
KZ012-1113N260	3.5	26	
KZ012-1113N280	3.5	28	
KZ012-1113N300	3.5	30	
KZ012-1113N320	3.5	32	
KZ012-1113N340	3.5	34	
KZ012-1113N360	3.5	36	
KZ012-1113N380	3.5	38	
KZ012-1113N400	3.5	40	
KZ012-1113N420	3.5	42	
KZ012-1113N440	3.5	44	
KZ012-1113N460	3.5	46	
KZ012-1113N480	3.5	48	
KZ012-1113N500	3.5	50	
KZ012-1113N520	3.5	52	
KZ012-1113N540	3.5	54	
KZ012-1113N560	3.5	56	
KZ012-1113N580	3.5	58	
KZ012-1113N600	3.5	60	
KZ012-1113N650	3.5	65	

螺钉系列Screw series

KZ012-1113N700	3.5	70	
KZ012-1113N750	3.5	75	
KZ012-1113N800	3.5	80	
KZ012-1113N850	3.5	85	
102.024.328	3.5	90	
备注: 钉帽形式为梅花型 Remark:The screw head type is Torx			

4.0锁定螺钉4.0 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1203N200	4.0	20	钛合金 Titanium Alloy
KZ012-1203N250	4.0	25	
KZ012-1203N300	4.0	30	
KZ012-1203N350	4.0	35	
KZ012-1203N400	4.0	40	
KZ012-1203N450	4.0	45	
KZ012-1203N500	4.0	50	
KZ012-1203N550	4.0	55	
KZ012-1203N600	4.0	60	
KZ012-1203N650	4.0	65	
KZ012-1203N700	4.0	70	
KZ012-1203N750	4.0	75	
KZ012-1203N800	4.0	80	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

5.0锁定螺钉5.0 Locking screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
102.025.300	5.0	14	钛合金 Titanium Alloy
102.025.301	5.0	16	
102.025.302	5.0	18	
KZ012-1123N200	5.0	20	
KZ012-1123N220	5.0	22	
KZ012-1123N240	5.0	24	
KZ012-1123N260	5.0	26	
KZ012-1123N280	5.0	28	
KZ012-1123N300	5.0	30	
KZ012-1123N320	5.0	32	
KZ012-1123N340	5.0	34	
KZ012-1123N360	5.0	36	
KZ012-1123N380	5.0	38	
KZ012-1123N400	5.0	40	
KZ012-1123N420	5.0	42	
KZ012-1123N440	5.0	44	
KZ012-1123N460	5.0	46	
KZ012-1123N480	5.0	48	
KZ012-1123N500	5.0	50	
KZ012-1123N550	5.0	55	
KZ012-1123N600	5.0	60	
KZ012-1123N650	5.0	65	
KZ012-1123N700	5.0	70	
KZ012-1123N750	5.0	75	
KZ012-1123N800	5.0	80	
KZ012-1123N850	5.0	85	
KZ012-1123N900	5.0	90	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列Screw series

6.5 锁定螺钉6.5 Locking Screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1203N500	6.5	50	钛合金 Titanium Alloy
KZ012-1203N550	6.5	55	
KZ012-1203N600	6.5	60	
KZ012-1203N650	6.5	65	
KZ012-1203N700	6.5	70	
KZ012-1203N750	6.5	75	
KZ012-1203N800	6.5	80	
KZ012-1203N850	6.5	85	
KZ012-1203N900	6.5	90	
KZ012-1203N950	6.5	95	
KZ012-1203N1000	6.5	100	
KZ012-1203N1050	6.5	105	
KZ012-1203N1100	6.5	110	
备注: 钉帽形式为梅花型 Remark:The screw head type is Torx			

7.3 锁定螺钉7.3 Locking Screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ012-1213N500	7.3	50	钛合金 Titanium Alloy
KZ012-1213N550	7.3	55	
KZ012-1213N600	7.3	60	
KZ012-1213N650	7.3	65	
KZ012-1213N700	7.3	70	
KZ012-1213N750	7.3	75	
KZ012-1213N800	7.3	80	
KZ012-1213N850	7.3	85	
KZ012-1213N900	7.3	90	
KZ012-1213N950	7.3	95	
KZ012-1213N1000	7.3	100	
KZ012-1213N1050	7.3	105	
KZ012-1213N1100	7.3	110	
备注: 钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列Screw series

6.5 锁定螺钉（空心）6.5 Locking Screw (hollow)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ012-1173N090	6.5	45	钛合金 Titanium Alloy
	KZ012-1173N100	6.5	50	
	KZ012-1173N110	6.5	55	
	KZ012-1173N120	6.5	60	
	KZ012-1173N130	6.5	65	
	KZ012-1173N140	6.5	70	
	KZ012-1173N150	6.5	75	
	KZ012-1173N160	6.5	80	
	KZ012-1173N170	6.5	85	
	KZ012-1173N180	6.5	90	
	KZ012-1173N190	6.5	95	
	KZ012-1173N200	6.5	100	
	KZ012-1173N210	6.5	105	
	KZ012-1173N220	6.5	110	

6.5 锁定螺钉（空心）6.5 Locking Screw (hollow)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ012-1183N131	6.5	65	钛合金 Titanium Alloy
	KZ012-1183N141	6.5	70	
	KZ012-1183N151	6.5	75	
	KZ012-1183N161	6.5	80	
	KZ012-1183N171	6.5	85	
	KZ012-1183N181	6.5	90	
	KZ012-1183N191	6.5	95	
	KZ012-1183N201	6.5	100	
	KZ012-1183N211	6.5	105	
	KZ012-1183N221	6.5	110	

1.5 非锁定接骨螺钉1.5 Cortical screws



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1013N060	1.5	6	钛合金 Titanium Alloy
KZ013-1013N070	1.5	7	
KZ013-1013N080	1.5	8	
KZ013-1013N090	1.5	9	
KZ013-1013N100	1.5	10	
KZ013-1013N110	1.5	11	
KZ013-1013N120	1.5	12	
KZ013-1013N130	1.5	13	
KZ013-1013N140	1.5	14	
KZ013-1013N150	1.5	15	
KZ013-1013N160	1.5	16	
KZ013-1013N180	1.5	18	
KZ013-1013N200	1.5	20	
KZ013-1013N220	1.5	22	
KZ013-1013N240	1.5	24	
备注: 钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列Screw series

2.0 非锁定接骨螺钉2.0 Cortical screws

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1033N060	2.0	6	钛合金 Titanium Alloy
	KZ013-1033N070	2.0	7	
	KZ013-1033N080	2.0	8	
	KZ013-1033N090	2.0	9	
	KZ013-1033N100	2.0	10	
	KZ013-1033N110	2.0	11	
	KZ013-1033N120	2.0	12	
	KZ013-1033N130	2.0	13	
	KZ013-1033N140	2.0	14	
	KZ013-1033N160	2.0	16	
	KZ013-1033N180	2.0	18	
	KZ013-1033N200	2.0	20	
	KZ013-1033N220	2.0	22	
	KZ013-1033N240	2.0	24	
	KZ013-1033N260	2.0	26	
	KZ013-1033N280	2.0	28	
	KZ013-1033N300	2.0	30	
	备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

2.5 非锁定接骨螺钉2.5 Cortical screws

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1053N060	2.5	6	钛合金 Titanium Alloy
	KZ013-1053N070	2.5	7	
	KZ013-1053N080	2.5	8	
	KZ013-1053N090	2.5	9	
	KZ013-1053N100	2.5	10	
	KZ013-1053N110	2.5	11	
	KZ013-1053N120	2.5	12	
	KZ013-1053N140	2.5	14	
	KZ013-1053N160	2.5	16	
	KZ013-1053N180	2.5	18	
	KZ013-1053N200	2.5	20	
	KZ013-1053N220	2.5	22	
	KZ013-1053N240	2.5	24	
	KZ013-1053N260	2.5	26	
	KZ013-1053N280	2.5	28	
	KZ013-1053N300	2.5	30	
	备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

2.7 非锁定接骨螺钉2.7 Cortical screws

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1063N060	2.7	6	
	KZ013-1063N070	2.7	7	
	KZ013-1063N080	2.7	8	
	KZ013-1063N090	2.7	9	
	KZ013-1063N100	2.7	10	
	KZ013-1063N110	2.7	11	
	KZ013-1063N120	2.7	12	

螺钉系列Screw series

	KZ013-1063N140	2.7	14	钛合金 Titanium Alloy
	KZ013-1063N160	2.7	16	
	KZ013-1063N180	2.7	18	
	KZ013-1063N200	2.7	20	
	KZ013-1063N220	2.7	22	
	KZ013-1063N240	2.7	24	
	KZ013-1063N260	2.7	26	
	KZ013-1063N280	2.7	28	
	KZ013-1063N300	2.7	30	
	KZ013-1063N320	2.7	32	
	KZ013-1063N340	2.7	34	
	KZ013-1063N360	2.7	36	
	KZ013-1063N380	2.7	38	
	KZ013-1063N400	2.7	40	
	备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

3.5 非锁定接骨螺钉3.5 Cortical screws



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1083N100	3.5	10	钛合金 Titanium Alloy
KZ013-1083N120	3.5	12	
KZ013-1083N140	3.5	14	
KZ013-1083N160	3.5	16	
KZ013-1083N180	3.5	18	
KZ013-1083N200	3.5	20	
KZ013-1083N220	3.5	22	
KZ013-1083N240	3.5	24	
KZ013-1083N260	3.5	26	
KZ013-1083N280	3.5	28	
KZ013-1083N300	3.5	30	
KZ013-1083N320	3.5	32	
KZ013-1083N340	3.5	34	
KZ013-1083N360	3.5	36	
KZ013-1083N380	3.5	38	
KZ013-1083N400	3.5	40	
KZ013-1083N420	3.5	42	
KZ013-1083N440	3.5	44	
KZ013-1083N450	3.5	45	
KZ013-1083N460	3.5	46	
KZ013-1083N480	3.5	48	
KZ013-1083N500	3.5	50	
KZ013-1083N520	3.5	52	
KZ013-1083N540	3.5	54	
KZ013-1083N550	3.5	55	
KZ013-1083N560	3.5	56	
KZ013-1083N580	3.5	58	
KZ013-1083N600	3.5	60	
KZ013-1083N650	3.5	65	
KZ013-1083N700	3.5	70	
KZ013-1083N750	3.5	75	
KZ013-1083N800	3.5	80	
KZ013-1083N850	3.5	85	
KZ013-1083N900	3.5	90	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列 Screw series

4.0 非锁定接骨螺钉 4.0 Cortical screws



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1163N100	4.0	10	钛合金 Titanium Alloy
KZ013-1163N120	4.0	12	
KZ013-1163N140	4.0	14	
KZ013-1163N160	4.0	16	
KZ013-1163N180	4.0	18	
KZ013-1163N200	4.0	20	
KZ013-1163N220	4.0	22	
KZ013-1163N240	4.0	24	
KZ013-1163N260	4.0	26	
KZ013-1163N280	4.0	28	
KZ013-1163N300	4.0	30	
KZ013-1163N320	4.0	32	
KZ013-1163N350	4.0	35	
KZ013-1163N360	4.0	36	
KZ013-1163N400	4.0	40	
KZ013-1163N450	4.0	45	
KZ013-1163N500	4.0	50	
KZ013-1163N550	4.0	55	
KZ013-1163N600	4.0	60	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

4.5 非锁定接骨螺钉 4.5 Cortical screws



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1103N200	4.5	20	钛合金 Titanium Alloy
KZ013-1103N220	4.5	22	
KZ013-1103N240	4.5	24	
KZ013-1103N260	4.5	26	
KZ013-1103N280	4.5	28	
KZ013-1103N300	4.5	30	
KZ013-1103N320	4.5	32	
KZ013-1103N340	4.5	34	
KZ013-1103N360	4.5	36	
KZ013-1103N380	4.5	38	
KZ013-1103N400	4.5	40	
KZ013-1103N420	4.5	42	
KZ013-1103N440	4.5	44	
KZ013-1103N460	4.5	46	
KZ013-1103N480	4.5	48	
KZ013-1103N500	4.5	50	
KZ013-1103N520	4.5	52	
KZ013-1103N540	4.5	54	
KZ013-1103N560	4.5	56	
KZ013-1103N580	4.5	58	
KZ013-1103N620	4.5	62	
KZ013-1103N640	4.5	64	
KZ013-1103N660	4.5	66	
KZ013-1103N680	4.5	68	
KZ013-1103N700	4.5	70	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

螺钉系列 Screw series

6.5 非锁定接骨螺钉 6.5 Cortical screws



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1193N250	6.5	25	钛合金 Titanium Alloy
KZ013-1193N300	6.5	30	
KZ013-1193N350	6.5	35	
KZ013-1193N400	6.5	40	
KZ013-1193N450	6.5	45	
KZ013-1193N500	6.5	50	
KZ013-1193N550	6.5	55	
KZ013-1193N600	6.5	60	
KZ013-1193N650	6.5	65	
KZ013-1193N700	6.5	70	
KZ013-1193N750	6.5	75	
KZ013-1193N800	6.5	80	
KZ013-1193N850	6.5	85	
KZ013-1193N900	6.5	90	
KZ013-1193N910	6.5	95	
KZ013-1193N920	6.5	100	
KZ013-1193N930	6.5	105	
KZ013-1193N940	6.5	110	
备注：钉帽形式为梅花型 Remark:The screw head type is Torx			

3.5 球头空心螺钉 3.5 Hollow ball head screw

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ031-1033N200	3.5	20	钛合金 Titanium Alloy
	KZ031-1033N220	3.5	22	
	KZ031-1033N240	3.5	24	
	KZ031-1033N260	3.5	26	
	KZ031-1033N280	3.5	28	
	KZ031-1033N300	3.5	30	
	KZ031-1033N320	3.5	32	
	KZ031-1033N340	3.5	34	
	KZ031-1033N360	3.5	36	
	KZ031-1033N380	3.5	38	
	KZ031-1033N400	3.5	40	
	KZ031-1033N420	3.5	42	
	KZ031-1033N440	3.5	44	
	KZ031-1033N460	3.5	46	
	KZ031-1033N480	3.5	48	
	KZ031-1033N500	3.5	50	

螺钉系列 Screw series

4.0 球头空心螺钉 4.0 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1093N200	4.0	20	钛合金 Titanium Alloy
KZ031-1093N220	4.0	22	
KZ031-1093N240	4.0	24	
KZ031-1093N260	4.0	26	
KZ031-1093N280	4.0	28	
KZ031-1093N300	4.0	30	
KZ031-1093N320	4.0	32	
KZ031-1093N340	4.0	34	
KZ031-1093N360	4.0	36	
KZ031-1093N380	4.0	38	
KZ031-1093N400	4.0	40	
KZ031-1093N420	4.0	42	
KZ031-1093N440	4.0	44	

4.0 球头空心螺钉 4.0 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1093N460	4.0	46	钛合金 Titanium Alloy
KZ031-1093N480	4.0	48	
KZ031-1093N500	4.0	50	
KZ031-1093N520	4.0	52	
KZ031-1093N540	4.0	54	
KZ031-1093N560	4.0	56	
KZ031-1093N580	4.0	58	
KZ031-1093N600	4.0	60	
KZ031-1093N620	4.0	62	
KZ031-1093N640	4.0	64	
KZ031-1093N660	4.0	66	
KZ031-1093N680	4.0	68	
KZ031-1093N700	4.0	70	

4.5 球头空心螺钉 4.5 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1043N200	4.5	20	钛合金 Titanium Alloy
KZ031-1043N220	4.5	22	
KZ031-1043N240	4.5	24	
KZ031-1043N260	4.5	26	
KZ031-1043N280	4.5	28	
KZ031-1043N300	4.5	30	
KZ031-1043N320	4.5	32	
KZ031-1043N340	4.5	34	
KZ031-1043N360	4.5	36	
KZ031-1043N380	4.5	38	
KZ031-1043N400	4.5	40	
KZ031-1043N420	4.5	42	
KZ031-1043N440	4.5	44	
KZ031-1043N460	4.5	46	
KZ031-1043N480	4.5	48	
KZ031-1043N500	4.5	50	
KZ031-1043N520	4.5	52	
KZ031-1043N540	4.5	54	
KZ031-1043N560	4.5	56	
KZ031-1043N580	4.5	58	
KZ031-1043N600	4.5	60	

螺钉系列 Screw series

6.5 球头空心螺钉 6.5 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1053H090	6.5	45	钛合金 Titanium Alloy
KZ031-1053H100	6.5	50	
KZ031-1053H110	6.5	55	
KZ031-1053H120	6.5	60	
KZ031-1053H130	6.5	65	
KZ031-1053H140	6.5	70	
KZ031-1053H150	6.5	75	
KZ031-1053H160	6.5	80	
KZ031-1053H170	6.5	85	
KZ031-1053H180	6.5	90	
KZ031-1053H190	6.5	95	
KZ031-1053H200	6.5	100	
KZ031-1053H210	6.5	105	
KZ031-1053H220	6.5	110	

6.5 球头空心螺钉 6.5 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1053S060	6.5	30	钛合金 Titanium Alloy
KZ031-1053S070	6.5	35	
KZ031-1053S080	6.5	40	
KZ031-1053S090	6.5	45	
KZ031-1053S100	6.5	50	
KZ031-1053S110	6.5	55	
KZ031-1053S120	6.5	60	
KZ031-1053S130	6.5	65	
KZ031-1053S140	6.5	70	
KZ031-1053S150	6.5	75	
KZ031-1053S160	6.5	80	
KZ031-1053S170	6.5	85	
KZ031-1053S180	6.5	90	
KZ031-1053S190	6.5	95	
KZ031-1053S200	6.5	100	
KZ031-1053S210	6.5	105	
KZ031-1053S220	6.5	110	

7.3 球头空心螺钉 7.3 Hollow ball head screw

Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1023S060	7.3	30	钛合金 Titanium Alloy
KZ031-1023S070	7.3	35	
KZ031-1023S080	7.3	40	
KZ031-1023S090	7.3	45	
KZ031-1023S100	7.3	50	
KZ031-1023S110	7.3	55	
KZ031-1023S120	7.3	60	
KZ031-1023S130	7.3	65	
KZ031-1023S140	7.3	70	
KZ031-1023S150	7.3	75	
KZ031-1023S160	7.3	80	
KZ031-1023S170	7.3	85	
KZ031-1023S180	7.3	90	
KZ031-1023S190	7.3	95	
KZ031-1023S200	7.3	100	
KZ031-1023S210	7.3	105	
KZ031-1023S220	7.3	110	

螺钉系列Screw series

7.3球头空心螺钉7.3 Hollow ball head screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ031-1023H090	7.3	45	钛合金 Titanium Alloy
KZ031-1023H100	7.3	50	
KZ031-1023H110	7.3	55	
KZ031-1023H120	7.3	60	
KZ031-1023H130	7.3	65	
KZ031-1023H140	7.3	70	
KZ031-1023H150	7.3	75	
KZ031-1023H160	7.3	80	
KZ031-1023H170	7.3	85	
KZ031-1023H180	7.3	90	
KZ031-1023H190	7.3	95	
KZ031-1023H200	7.3	100	
KZ031-1023H210	7.3	105	
KZ031-1023H220	7.3	110	

2.5 无头加压空心螺钉II2.5 Headless compression hollow screw II



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ033-1013N080	2.5	8	钛合金 Titanium Alloy
KZ033-1013N100	2.5	10	
KZ033-1013N120	2.5	12	
KZ033-1013N140	2.5	14	
KZ033-1013N160	2.5	16	
KZ033-1013N180	2.5	18	
KZ033-1013N200	2.5	20	


3.0 无头加压空心螺钉II3.0 Headless compression hollow screw II



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ033-1023N080	3.0	8	钛合金 Titanium Alloy
KZ033-1023N100	3.0	10	
KZ033-1023N120	3.0	12	
KZ033-1023N140	3.0	14	
KZ033-1023N160	3.0	16	
KZ033-1023N180	3.0	18	
KZ033-1023N200	3.0	20	
KZ033-1023N220	3.0	22	
KZ033-1023N240	3.0	24	
KZ033-1023N260	3.0	26	
KZ033-1023N280	3.0	28	
KZ033-1023N300	3.0	30	

螺钉系列Screw series

3.5 无头加压空心螺钉II3.5 Headless compression hollow screw II



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ033-1033N080	3.5	8	钛合金 Titanium Alloy
KZ033-1033N100	3.5	10	
KZ033-1033N120	3.5	12	
KZ033-1033N140	3.5	14	
KZ033-1033N160	3.5	16	
KZ033-1033N180	3.5	18	
KZ033-1033N200	3.5	20	
KZ033-1033N220	3.5	22	
KZ033-1033N240	3.5	24	
KZ033-1033N260	3.5	26	
KZ033-1033N280	3.5	28	
KZ033-1033N300	3.5	30	

4.0 无头加压空心螺钉II4.0 Headless compression hollow screw II



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ033-1043N160	4.0	16	钛合金 Titanium Alloy
KZ033-1043N180	4.0	18	
KZ033-1043N200	4.0	20	
KZ033-1043N220	4.0	22	
KZ033-1043N240	4.0	24	
KZ033-1043N260	4.0	26	
KZ033-1043N280	4.0	28	
KZ033-1043N300	4.0	30	
KZ033-1043N320	4.0	32	
KZ033-1043N340	4.0	34	
KZ033-1043N360	4.0	36	
KZ033-1043N380	4.0	38	
KZ033-1043N400	4.0	40	
KZ033-1043N420	4.0	42	
KZ033-1043N440	4.0	44	
KZ033-1043N460	4.0	46	
KZ033-1043N480	4.0	48	
KZ033-1043N500	4.0	50	

5.0 无头加压空心螺钉II5.0 Headless compression hollow screw II



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ033-1053N200	5.0	20	钛合金 Titanium Alloy
KZ033-1053N250	5.0	25	
KZ033-1053N300	5.0	30	
KZ033-1053N350	5.0	35	
KZ033-1053N400	5.0	40	
KZ033-1053N450	5.0	45	
KZ033-1053N500	5.0	50	

螺钉系列Screw series

6.5 无头加压空心螺钉II6.5 Headless compression hollow screw II

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ033-1063N450	6.5	45	钛合金 Titanium Alloy
	KZ033-1063N500	6.5	50	
	KZ033-1063N550	6.5	55	
	KZ033-1063N600	6.5	60	
	KZ033-1063N650	6.5	65	
	KZ033-1063N700	6.5	70	
	KZ033-1063N750	6.5	75	

6.5 锥形空心螺钉II6.5 Conical Hollow Screw II

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ034-1013S070	6.5	35	钛合金 Titanium Alloy
	KZ034-1013S080	6.5	40	
	KZ034-1013S090	6.5	45	
	KZ034-1013S100	6.5	50	
	KZ034-1013S110	6.5	55	
	KZ034-1013S120	6.5	60	
	KZ034-1013S130	6.5	65	
	KZ034-1013S140	6.5	70	
	KZ034-1013S150	6.5	75	
	KZ034-1013S160	6.5	80	
	KZ034-1013S170	6.5	85	
	KZ034-1013S180	6.5	90	
	KZ034-1013S190	6.5	95	
	KZ034-1013S200	6.5	100	
	KZ034-1013S210	6.5	105	
	KZ034-1013S220	6.5	110	

6.5 锥形空心螺钉II6.5 Conical Hollow Screw II

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ034-1013H090	6.5	45	钛合金 Titanium Alloy
	KZ034-1013H100	6.5	50	
	KZ034-1013H110	6.5	55	
	KZ034-1013H120	6.5	60	
	KZ034-1013H130	6.5	65	
	KZ034-1013H140	6.5	70	
	KZ034-1013H150	6.5	75	
	KZ034-1013H160	6.5	80	
	KZ034-1013H170	6.5	85	
	KZ034-1013H180	6.5	90	
	KZ034-1013H190	6.5	95	
	KZ034-1013H200	6.5	100	
	KZ034-1013H210	6.5	105	
	KZ034-1013H220	6.5	110	

螺钉系列Screw series

3.5 金属接骨螺钉（HB半牙）3.5 Metal bone screw (HB half-tooth)



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1323H120	3.5	12	钛合金 Titanium Alloy
KZ013-1323H140	3.5	14	
KZ013-1323H160	3.5	16	
KZ013-1323H180	3.5	18	
KZ013-1323H200	3.5	20	
KZ013-1323H220	3.5	22	
KZ013-1323H240	3.5	24	
KZ013-1323H260	3.5	26	
KZ013-1323H280	3.5	28	
KZ013-1323H300	3.5	30	
KZ013-1323H320	3.5	32	
KZ013-1323H340	3.5	34	
KZ013-1323H360	3.5	36	
KZ013-1323H380	3.5	38	
KZ013-1323H400	3.5	40	
KZ013-1323H420	3.5	42	
KZ013-1323H440	3.5	44	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal			

4.0 金属接骨螺钉（HB半牙）4.0 Metal bone screw (HB half-tooth)



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ013-1333H120	4.0	12	钛合金 Titanium Alloy
KZ013-1333H140	4.0	14	
KZ013-1333H160	4.0	16	
KZ013-1333H180	4.0	18	
KZ013-1333H200	4.0	20	
KZ013-1333H220	4.0	22	
KZ013-1333H240	4.0	24	
KZ013-1333H260	4.0	26	
KZ013-1333H280	4.0	28	
KZ013-1333H300	4.0	30	
KZ013-1333H320	4.0	32	
KZ013-1333H340	4.0	34	
KZ013-1333H360	4.0	36	
KZ013-1333H380	4.0	38	
KZ013-1333H400	4.0	40	
KZ013-1333H420	4.0	42	
KZ013-1333H440	4.0	44	
KZ013-1333H460	4.0	46	
KZ013-1333H480	4.0	48	
KZ013-1333H500	4.0	50	
KZ013-1333H520	4.0	52	
KZ013-1333H540	4.0	54	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal			

螺钉系列Screw series

4.5 金属接骨螺钉（HB半牙）4.5 Metal bone screw (HB half-tooth)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1343H200	4.5	20	钛合金 Titanium Alloy
	KZ013-1343H220	4.5	22	
	KZ013-1343H240	4.5	24	
	KZ013-1343H260	4.5	26	
	KZ013-1343H280	4.5	28	
	KZ013-1343H300	4.5	30	
	KZ013-1343H320	4.5	32	
	KZ013-1343H340	4.5	34	
	KZ013-1343H360	4.5	36	
	KZ013-1343H380	4.5	38	
	KZ013-1343H400	4.5	40	
	KZ013-1343H420	4.5	42	
	KZ013-1343H440	4.5	44	
	KZ013-1343H460	4.5	46	
	KZ013-1343H480	4.5	48	
	KZ013-1343H500	4.5	50	
	KZ013-1343H520	4.5	52	
	KZ013-1343H540	4.5	54	
	KZ013-1343H560	4.5	56	
	KZ013-1343H580	4.5	58	
	KZ013-1343H600	4.5	60	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal				

6.5 金属接骨螺钉（HB半牙）6.5 Metal bone screw (HB half-tooth)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1443H350	6.5	35	钛合金 Titanium Alloy
	KZ013-1443H400	6.5	40	
	KZ013-1443H450	6.5	45	
	KZ013-1443H500	6.5	50	
	KZ013-1443H550	6.5	55	
	KZ013-1443H600	6.5	60	
	KZ013-1443H650	6.5	65	
	KZ013-1443H700	6.5	70	
	KZ013-1443H750	6.5	75	
	KZ013-1443H800	6.5	80	
	KZ013-1443H850	6.5	85	
	KZ013-1443H900	6.5	90	
	KZ013-1443H950	6.5	95	
	KZ013-1443H1000	6.5	100	
	KZ013-1443H1050	6.5	105	
	KZ013-1443H1100	6.5	110	
	KZ013-1443H1150	6.5	115	
	KZ013-1443H1200	6.5	120	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal				

螺钉系列Screw series

4.0 金属接骨螺钉（HB全牙）4.0 Metal bone screws (HB full teeth)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1353N120	4.0	12	钛合金 Titanium Alloy
	KZ013-1353N140	4.0	14	
	KZ013-1353N160	4.0	16	
	KZ013-1353N180	4.0	18	
	KZ013-1353N200	4.0	20	
	KZ013-1353N220	4.0	22	
	KZ013-1353N240	4.0	24	
	KZ013-1353N260	4.0	26	
	KZ013-1353N280	4.0	28	
	KZ013-1353N300	4.0	30	
	KZ013-1353N320	4.0	32	
	KZ013-1353N340	4.0	34	
	KZ013-1353N360	4.0	36	
	KZ013-1353N380	4.0	38	
	KZ013-1353N400	4.0	40	
	KZ013-1353N420	4.0	42	
	KZ013-1353N440	4.0	44	
	KZ013-1353N460	4.0	46	
	KZ013-1353N480	4.0	48	
	KZ013-1353N500	4.0	50	
	KZ013-1353N520	4.0	52	
	KZ013-1353N540	4.0	54	
	KZ013-1353N560	4.0	56	
	KZ013-1353N580	4.0	58	
	KZ013-1353N600	4.0	60	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal				

6.5 金属接骨螺钉（HB全牙）6.5 Metal bone screws (HB full teeth)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1363N300	6.5	30	钛合金 Titanium Alloy
	KZ013-1363N350	6.5	35	
	KZ013-1363N400	6.5	40	
	KZ013-1363N450	6.5	45	
	KZ013-1363N500	6.5	50	
	KZ013-1363N550	6.5	55	
	KZ013-1363N600	6.5	60	
	KZ013-1363N650	6.5	65	
	KZ013-1363N700	6.5	70	
	KZ013-1363N750	6.5	75	
	KZ013-1363N800	6.5	80	
	KZ013-1363N850	6.5	85	
	KZ013-1363N900	6.5	90	
	KZ013-1363N950	6.5	95	
	KZ013-1363N1000	6.5	100	
	KZ013-1363N1050	6.5	105	
	KZ013-1363N1100	6.5	110	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal				

螺钉系列Screw series

2.0 金属接骨螺钉（HA全牙）2.0 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1373N060	2.0	6	钛合金 Titanium Alloy
	KZ013-1373N080	2.0	8	
	KZ013-1373N100	2.0	10	
	KZ013-1373N120	2.0	12	
	KZ013-1373N140	2.0	14	
	KZ013-1373N160	2.0	16	
	KZ013-1373N180	2.0	18	
	KZ013-1373N200	2.0	20	
	KZ013-1373N220	2.0	22	
	备注：钉帽形式为梅花型和四方型 Remark:The screw head types are Torx and square			

2.5 金属接骨螺钉（HA全牙）2.5 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1383N060	2.5	6	钛合金 Titanium Alloy
	KZ013-1383N080	2.5	8	
	KZ013-1383N100	2.5	10	
	KZ013-1383N120	2.5	12	
	KZ013-1383N140	2.5	14	
	KZ013-1383N160	2.5	16	
	KZ013-1383N180	2.5	18	
	KZ013-1383N200	2.5	20	
	KZ013-1383N220	2.5	22	
	备注：钉帽形式为梅花型和四方型 Remark:The screw head types are Torx and square			

2.7 金属接骨螺钉（HA全牙）2.7 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1393N060	2.7	6	钛合金 Titanium Alloy
	KZ013-1393N080	2.7	8	
	KZ013-1393N100	2.7	10	
	KZ013-1393N120	2.7	12	
	KZ013-1393N140	2.7	14	
	KZ013-1393N160	2.7	16	
	KZ013-1393N180	2.7	18	
	KZ013-1393N200	2.7	20	
	KZ013-1393N220	2.7	22	
	KZ013-1393N240	2.7	24	
	KZ013-1393N260	2.7	26	
	KZ013-1393N280	2.7	28	
	KZ013-1393N300	2.7	30	
	备注：钉帽形式为梅花型和四方型 Remark:The screw head types are Torx and square			

螺钉系列Screw series

3.5 金属接骨螺钉（HA全牙）3.5 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1403N120	3.5	12	钛合金 Titanium Alloy
	KZ013-1403N140	3.5	14	
	KZ013-1403N160	3.5	16	
	KZ013-1403N180	3.5	18	
	KZ013-1403N200	3.5	20	
	KZ013-1403N220	3.5	22	
	KZ013-1403N240	3.5	24	
	KZ013-1403N260	3.5	26	
	KZ013-1403N280	3.5	28	
	KZ013-1403N300	3.5	30	
	KZ013-1403N320	3.5	32	
	KZ013-1403N340	3.5	34	
	KZ013-1403N360	3.5	36	
	KZ013-1403N380	3.5	38	
	KZ013-1403N400	3.5	40	
	KZ013-1403N420	3.5	42	
	KZ013-1403N440	3.5	44	
	备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal			

4.0 金属接骨螺钉（HA全牙）4.0 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1413N140	4.0	14	钛合金 Titanium Alloy
	KZ013-1413N160	4.0	16	
	KZ013-1413N180	4.0	18	
	KZ013-1413N200	4.0	20	
	KZ013-1413N220	4.0	22	
	KZ013-1413N240	4.0	24	
	KZ013-1413N260	4.0	26	
	KZ013-1413N280	4.0	28	
	KZ013-1413N300	4.0	30	
	KZ013-1413N320	4.0	32	
	KZ013-1413N340	4.0	34	
	KZ013-1413N360	4.0	36	
	KZ013-1413N380	4.0	38	
	KZ013-1413N400	4.0	40	
	KZ013-1413N420	4.0	42	
	KZ013-1413N440	4.0	44	
	KZ013-1413N460	4.0	46	
	KZ013-1413N480	4.0	48	
	KZ013-1413N500	4.0	50	
	KZ013-1413N520	4.0	52	
	KZ013-1413N540	4.0	54	
	备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal			

螺钉系列Screw series

4.5 金属接骨螺钉（HA全牙）4.5 Metal bone screws (HA full dental)

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ013-1453N180	4.5	18	钛合金 Titanium Alloy
	KZ013-1453N200	4.5	20	
	KZ013-1453N220	4.5	22	
	KZ013-1453N240	4.5	24	
	KZ013-1453N260	4.5	26	
	KZ013-1453N280	4.5	28	
	KZ013-1453N300	4.5	30	
	KZ013-1453N320	4.5	32	
	KZ013-1453N340	4.5	34	
	KZ013-1453N360	4.5	36	
	KZ013-1453N380	4.5	38	
	KZ013-1453N400	4.5	40	
	KZ013-1453N420	4.5	42	
	KZ013-1453N440	4.5	44	
	KZ013-1453N460	4.5	46	
	KZ013-1453N480	4.5	48	
	KZ013-1453N500	4.5	50	
	KZ013-1453N520	4.5	52	
	KZ013-1453N540	4.5	54	
	KZ013-1453N560	4.5	56	
	KZ013-1453N580	4.5	58	
	KZ013-1453N600	4.5	60	
备注：钉帽形式为梅花型和六角型 Remark:The screw head types are Torx and hexagonal				

3.0 双头加压空心螺钉3.0 Double-ended pressure hollow screw

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ033-1103N140	3.0	14	钛合金 Titanium Alloy
	KZ033-1103N160	3.0	16	
	KZ033-1103N180	3.0	18	
	KZ033-1103N200	3.0	20	
	KZ033-1103N220	3.0	22	
	KZ033-1103N240	3.0	24	
	KZ033-1103N260	3.0	26	
	KZ033-1103N280	3.0	28	
	KZ033-1103N300	3.0	30	

螺钉系列Screw series

4.0 双头加压空心螺钉4.0 Double-ended pressure hollow screw

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ033-1113N140	4.0	14	钛合金 Titanium Alloy
	KZ033-1113N140	4.0	16	
	KZ033-1113N140	4.0	18	
	KZ033-1113N140	4.0	20	
	KZ033-1113N140	4.0	22	
	KZ033-1113N140	4.0	24	
	KZ033-1113N140	4.0	26	
	KZ033-1113N140	4.0	28	
	KZ033-1113N140	4.0	30	

3.5 空心接骨螺钉3.5 Hollow bone screw

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ034-1203N180	3.5	18	钛合金 Titanium Alloy
	KZ034-1203N200	3.5	20	
	KZ034-1203N220	3.5	22	
	KZ034-1203N240	3.5	24	
	KZ034-1203N260	3.5	26	
	KZ034-1203N280	3.5	28	
	KZ034-1203N300	3.5	30	
	KZ034-1203N320	3.5	32	
	KZ034-1203N340	3.5	34	
	KZ034-1203N360	3.5	36	
	KZ034-1203N380	3.5	38	
	KZ034-1203N400	3.5	40	
	KZ034-1203N420	3.5	42	
	KZ034-1203N440	3.5	44	
	KZ034-1203N460	3.5	46	
	KZ034-1203N480	3.5	48	
	KZ034-1203N500	3.5	50	
	KZ034-1203N520	3.5	52	
	KZ034-1203N540	3.5	54	
	KZ034-1203N560	3.5	56	
	KZ034-1203N580	3.5	58	
	KZ034-1203N560	3.5	60	

螺钉系列Screw series

4.0 空心接骨螺钉4.0 Hollow bone screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ034-1213N200	4.0	20	钛合金 Titanium Alloy
KZ034-1213N220	4.0	22	
KZ034-1213N240	4.0	24	
KZ034-1213N260	4.0	26	
KZ034-1213N280	4.0	28	
KZ034-1213N300	4.0	30	
KZ034-1213N320	4.0	32	
KZ034-1213N340	4.0	34	
KZ034-1213N360	4.0	36	
KZ034-1213N380	4.0	38	
KZ034-1213N400	4.0	40	
KZ034-1213N420	4.0	42	
KZ034-1213N440	4.0	44	
KZ034-1213N460	4.0	46	
KZ034-1213N480	4.0	48	
KZ034-1213N500	4.0	50	
KZ034-1213N520	4.0	52	
KZ034-1213N540	4.0	54	
KZ034-1213N560	4.0	56	
KZ034-1213N580	4.0	58	
KZ034-1213N560	4.0	60	


4.5 空心接骨螺钉4.5 Hollow bone screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ034-1223N300	4.5	30	钛合金 Titanium Alloy
KZ034-1223N320	4.5	32	
KZ034-1223N340	4.5	34	
KZ034-1223N360	4.5	36	
KZ034-1223N380	4.5	38	
KZ034-1223N400	4.5	40	
KZ034-1223N420	4.5	42	
KZ034-1223N440	4.5	44	
KZ034-1223N460	4.5	46	
KZ034-1223N480	4.5	48	
KZ034-1223N500	4.5	50	
KZ034-1223N520	4.5	52	
KZ034-1223N540	4.5	54	
KZ034-1223N560	4.5	56	
KZ034-1223N580	4.5	58	
KZ034-1223N560	4.5	60	
KZ034-1223N562	4.5	62	
KZ034-1223N564	4.5	64	

螺钉系列Screw series

6.5 空心接骨螺钉6.5 Hollow bone screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ034-1233N500	6.5	50	钛合金 Titanium Alloy
KZ034-1233N550	6.5	55	
KZ034-1233N600	6.5	60	
KZ034-1233N650	6.5	65	
KZ034-1233N700	6.5	70	
KZ034-1233N750	6.5	75	
KZ034-1233N800	6.5	80	
KZ034-1233N850	6.5	85	
KZ034-1233N900	6.5	90	
KZ034-1233N950	6.5	95	
KZ034-1233N1000	6.5	100	
KZ034-1233N1050	6.5	105	
KZ034-1233N1100	6.5	110	
KZ034-1233N1150	6.5	115	
KZ034-1233N1200	6.5	120	

7.3 空心接骨螺钉7.3 Hollow bone screw



Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
KZ034-1243N500	7.5	50	钛合金 Titanium Alloy
KZ034-1243N550	7.5	55	
KZ034-1243N600	7.5	60	
KZ034-1243N650	7.5	65	
KZ034-1243N700	7.5	70	
KZ034-1243N750	7.5	75	
KZ034-1243N800	7.5	80	
KZ034-1243N850	7.5	85	
KZ034-1243N900	7.5	90	
KZ034-1243N950	7.5	95	
KZ034-1243N1000	7.5	100	
KZ034-1243N1050	7.5	105	
KZ034-1243N1100	7.5	110	
KZ034-1243N1150	7.5	115	
KZ034-1243N1200	7.5	120	

螺钉系列 Screw series

空心双头加压螺钉 Hollow double-ended pressure screw

	Product No. 产品编号	Thread Diameter 螺纹直径	Length 长度 (mm)	Material 材质
	KZ032-1013N050	6.0	50	钛合金 Titanium Alloy
	KZ032-1013N055	6.0	55	
	KZ032-1013N060	6.0	60	
	KZ032-1013N065	6.0	65	
	KZ032-1013N070	6.0	70	
	KZ032-1013N075	6.0	75	
	KZ032-1013N080	6.0	80	
	KZ032-1013N085	6.0	85	
	KZ032-1013N090	6.0	90	
	KZ032-1013N095	6.0	95	
	KZ032-1013N100	6.0	100	
	KZ032-1013N105	6.0	105	
	KZ032-1013N110	6.0	110	

金属缆索系统 Metal cable system

1.0 无针型缆索 I 1.0 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1103N400	1.0	400	钛合金 Titanium Alloy
KZ041-1103N450	1.0	450	
KZ041-1103N500	1.0	500	
KZ041-1103N550	1.0	550	
KZ041-1103N600	1.0	600	
KZ041-1103N650	1.0	650	
KZ041-1103N700	1.0	700	

1.2 无针型缆索 I 1.2 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1123N400	1.2	400	钛合金 Titanium Alloy
KZ041-1123N450	1.2	450	
KZ041-1123N500	1.2	500	
KZ041-1123N550	1.2	550	
KZ041-1123N600	1.2	600	
KZ041-1123N650	1.2	650	
KZ041-1123N700	1.2	700	

1.4 无针型缆索 I 1.4 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1143N400	1.4	400	钛合金 Titanium Alloy
KZ041-1143N450	1.4	450	
KZ041-1143N500	1.4	500	
KZ041-1143N550	1.4	550	
KZ041-1143N600	1.4	600	
KZ041-1143N650	1.4	650	
KZ041-1143N700	1.4	700	
KZ041-1143N800	1.4	800	

1.6 无针型缆索 I 1.6 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1163N400	1.6	400	钛合金 Titanium Alloy
KZ041-1163N450	1.6	450	
KZ041-1163N500	1.6	500	
KZ041-1163N550	1.6	550	
KZ041-1163N600	1.6	600	
KZ041-1163N650	1.6	650	
KZ041-1163N700	1.6	700	

金属缆索系统

Metal cable system

1.8 无针型缆索 I

1.8 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1183N400	1.8	400	钛合金 Titanium Alloy
KZ041-1183N450	1.8	450	
KZ041-1183N500	1.8	500	
KZ041-1183N550	1.8	550	
KZ041-1183N600	1.8	600	
KZ041-1183N650	1.8	650	
KZ041-1183N700	1.8	700	
KZ041-1183N800	1.8	800	

2.0 无针型缆索 I

2.0 Needleless Cable I



Product No. 产品编号	Diameter 直径	Length 长度(mm)	Material 材质
KZ041-1203N400	2	400	钛合金 Titanium Alloy
KZ041-1203N450	2	450	
KZ041-1203N500	2	500	
KZ041-1203N550	2	550	
KZ041-1203N600	2	600	
KZ041-1203N650	2	650	
KZ041-1203N700	2	700	

压紧式锁紧扣

Compression locking buckle



Product No. 产品编号	Diameter 直径	Material 材质
KZ042-1113N000	1.1	钛合金 Titanium Alloy
KZ042-1133N000	1.3	
KZ042-1153N000	1.5	
KZ042-1173N000	1.7	
KZ042-1193N000	1.9	
KZ042-1213N000	2.1	